

**FORCEPS, HEMOSTATS, CLAMPS****PG131114A : LANDRENEAU PENNINGTON GRASPING FORCEPS**

Straight, 10.2 mm tip, 5-1/8" (13 cm)

**PG181220 : ADSON DRESSING FORCEPS**

Serrated, 4 3/4" (12 cm)

**PG181223 : ADSON TISSUE FORCEPS**

1 x 2 Teeth, 4 3/4" (12 cm)

**PG181225 : ADSON TISSUE FORCEPS**

2 x 3 Teeth, 5" (13 cm)

**PG181228 : ADSON DRESSING FORCEPS**

Serrated, 1 x 2 Teeth, 4 3/4" (12 cm)

**PG181235 : BROWN-ADSON TISSUE FORCEPS**

9 x 9 Teeth, 4 3/4" (12 cm)



**PG181237 : EWALD DRESSING FORCEPS**

Serrated, 4 3/4" (12 cm)



**PG181240 : EWALD TISSUE FORCEPS**

1 x 2 Teeth, 4 3/4" (12 cm)



**PG181255 : BROWN TISSUE FORCEPS**

Straight, 9 x 9 Teeth, 6" (16 cm)



**PG181360 : RUSSIAN TISSUE FORCEPS**

6" (15 cm)



**PG181365 : RUSSIAN TISSUE FORCEPS**

8" (21 cm)



**PG181370 : RUSSIAN TISSUE FORCEPS**

10" (26 cm)




**PG181410 : CUSHING DRESSING FORCEPS**

Serrated, 7" (18 cm)



**PG183010 : BABCOCK TISSUE FORCEPS**

6 1/4" (16 cm)



**PG183020 : BABCOCK TISSUE FORCEPS**

7 7/8" (20 cm)



**PG211116A : LANDRENEAU MASHER GRASPING FORCEPS**

Straight, 7.8 mm tip, 8-1/2" (22 cm)



**PG211117A : LANDRENEAU MASHER GRASPING FORCEPS**

Lightly curved, 7.8 mm tip, 8-1/2" (22 cm)



**PG211118A : LANDRENEAU MASHER GRASPING FORCEPS**

Straight, 13.8 mm tip, 8-1/2" (22 cm)



**PG211245A : LANDRENEAU LAWRENCE FINE DISSECTING FORCEPS**

45 degree curve, locking, 21 cm shaft length, 15-3/8" (39 cm) overall length. Not designed for use with a 10 mm cannula.



**PG211290A : LANDRENEAU MIXTER DISSECTING FORCEPS**

Curved, 2.7 mm tip, 8-1/2" (22 cm)




**PG211603A : LANDRENEAU DEBAKEY GRASPING FORCEPS**

Straight, 7 mm tip, 8-1/2" (22 cm)



**PG342557 : KANTROWITZ FORCEPS**

Right Angle, Serrated, 3/4" (1.75cm) tip length, 2mm tip width, 7 7/8" (20cm)




**PG342558 : KANTROWITZ RIGHT ANGLE DISSECTING FORCEPS**

Serrated, 9 1/2" (25 cm)



**PG342559 : KANTROWITZ FORCEPS**

Right Angle, Serrated, 3/4" (1.95cm) tip length, 2mm tip width, 10 5/8" (27cm)



**PG342560 : DEBAKEY FORCEPS**

C Curved, 3/8" (1cm) tip length, 8 1/8" (20.5cm)

