


FORCEPS, HEMOSTATS, CLAMPS**PG342561 : DEBAKEY FORCEPS**


D Curved, 1/2" (1.4cm) tip length, 2.8mm tip width, 8 1/4" (21cm)

**PG342562 : DEBAKEY FORCEPS**

C Curved, 5/8" (1.5cm) tip length, 2.8mm tip width, 9 5/8" (24.5cm)

**PG342563 : DEBAKEY FORCEPS**


D Curved, 7/8" (2.2cm) tip length, 9 1/2" (24cm)

**PG342565 : DEBAKEY FORCEPS**


D Curved, 3/4" (2cm) tip length, 10 1/4" (26cm)

**PG342566 : DEBAKEY THORACIC DISSECTING FORCEPS**

Angled, 7 5/8" (20 cm)


**PG342567 : DEBAKEY FORCEPS**

Right Angle, Longitudinal Serrations, Cross-serrated tip, 7/8" (2.3cm) tip length, 9" (23cm)




PG342568 : DEBAKEY FORCEPS

Right Angle, Longitudinal Serrations, Cross-serrated tip, 7/8" (2.14cm) tip length, 10 1/4" (26cm)




PG342569 : DEBAKEY FORCEPS

Right Angle, Longitudinal Serrations, Cross-serrated tip, 7/8" (2.15cm) tip length, 11 3/8" (29cm)




PG342571 : SAROT THORACIC ARTERY FORCEPS

Curved, Serrated, 1/2" (1.3cm) tip length, 1.6mm tip width, 7 5/8" (19.5cm)




PG342572 : SAROT THORACIC ARTERY FORCEPS

Curved, Serrated, 2 1/4" (5.7cm) tip length, 2mm tip width, 9" (23cm)




PG342575 : ROOS SCALENECTOMY FORCEPS

Right Angle, Longitudinal Serrations, 1/2" (1.2cm) tip length, 1.5mm tip width, 11" (28cm)



PG342580 : WILLAUER FORCEPS

Strongly Curved, Longitudinal Serrations, 2" (5cm) tip length, 9 1/2" (24cm)




PG342583 : SAWTELL FORCEPS

Full Curved, Serrated, 1 1/8" (2.7cm) tip length, 1.7mm tip width, 7 1/4" (18.5cm)




PG342585 : PILLING FORCEPS

Slight Curved, Serrated, 2" (5.1cm) tip length, 2mm tip width, 9 1/4" (23.5cm)




PG342587 : BENGOLEA FORCEPS

Slight Curved, Serrated, 2mm tip width, 8 1/2" (21.5cm)



PG342590 : VANDERBILT VESSEL FORCEPS

Slight Curved, Serrated, 2" (5cm) tip length, 2mm tip width, 8 5/8" (22cm)




PG342594 : LEJEUNE FORCEPS

Slight Curved, Serrated, 1 3/8" (3.5cm) tip length, 1.8mm tip width, 9 3/8" (24cm)



PG342596 : GEMINI FORCEPS

Curved, Serrated, 1/2" (1.23cm) tip length, 1.3mm tip width, 7 1/4" (18.5cm)



PG342597 : GEMINI FORCEPS

Curved, Serrated, 3/4" (1.9cm) tip length, 1.8mm tip width, 9" (23cm)




PG342598 : GEMINI THORACIC DISSECTING FORCEPS

Curved, Serrated, 11 1/8" (28 cm)




PG342600 : MIXER THORACIC FORCEPS

Curved, DeBakey Jaw, 1" (2.4cm) tip length, 2.8mm tip width, 9 1/4" (23.5cm)




PG342601 : MIXER FORCEPS

Strongly Curved, Serrated, 3/4" (2cm) tip length, 8 3/4" (22cm)




PG342603 : MIXER FORCEPS

Half Curved, Serrated, 5/8" (1.65cm) tip length, 2.3mm tip width, 9 1/8" (23cm)



PG342605 : MIXER FORCEPS

Strongly Curved, DeBakey Jaw, 1 3/8" (3.5cm) tip length, 11" (28cm)



PG342606 : MIXER FORCEPS

Strongly Curved, Serrated (3cm) tip length, 10 3/4" (27.5cm)

