

GENERAL**PG162415 : FERGUSON-FRAZIER SUCTION TUBE # 5**

12 Fr. tip, 4-7/8" (12.4 cm)

**PG162416 : FERGUSON-FRAZIER SUCTION TUBE # 6**


14 Fr. tip, 4-7/8" (12.4 cm)

**PG162421 : FERGUSON-FRAZIER SUCTION TUBE**

Ferguson Stylet # 0

**PG163000 : WISHBONE RETRACTOR SET**

Set contains 163100, 163102, 163103, 163104 and 163106

**PG163100 : WISHBONE PILLING TABLE CLAMP****PG163104 : WISHBONE SPLIT BAR**

15 X 17 inches



PG163108 : WISHBONE CLOSED UNIVERSAL CLAMP



PG163109 : WISHBONE OPEN UNIVERSAL CLAMP



PG163200 : WISHBONE SPLANCHNIC RETRACTOR

Curved, 2" (5cm) x 4 1/4" (11.4cm)



PG163201 : WISHBONE SPLANCHNIC RETRACTOR

Curved, 2" (5cm) x 6 1/4" (16.0cm)



PG163203 : WISHBONE RENAL VEIN RETRACTOR

Curved, 1 x 6, 7" (18cm)



PG163211 : WISHBONE MAYO BLADE STANDARD

Curved, 2 5/8" (6.7cm) x 2.25



PG163212 : WISHBONE MAYO BLADE WIDE

Curved, 3 1/2" (8.9cm) x 2.75

PG163213 : WISHBONE RICHARDSON SWIVEL BLADE

Curved, 1" (2.5cm) x 1.75



PG163400 : WISHBONE BLADE HOLDER SWIVEL HANDLE

6 1/4" (16cm)



PG163401 : WISHBONE BLADE HOLDER SWIVEL KNOB

9 1/4" (23.5cm)



PG163406 : WISHBONE FENCE SWIVEL BLADE

Curved, 4" (10.0cm), 6



PG163422 : WISHBONE MAYO WIDE SWIVEL BLADE

Curved, 3 1/2" (9.0cm), 2.75




PG163424 : WISHBONE DEAYER SWIVEL BLADE

Curved, 2 3/4" (7cm), 4.5



PG163436 : WISHBONE LECHNER-HARRINGTON LIVER BLADE

1 5/8" (4.0cm), 5, 10 1/4" (26cm)





PG163701 : RETRACTOR FIBEROPTIC LIGHT BEAM



PG164110 : DEAVER RETRACTOR

3/4" (1.9cm) width blade, 6 7/8" (17.5cm)



PG164112 : DEAVER RETRACTOR

7/8" (2.2cm) width blade, 7 7/8" (20.0cm)



PG164115 : DEAVER RETRACTOR

1" (2.5cm) width blade, 8 7/8" (22.5cm)



PG164118 : DEAVER RETRACTOR, #1

1" (2.5cm) width blade, 8 1/8" (20.5cm)

