

HORIZON™ METAL LIGATION SYSTEM - CLIPS**WK533837 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small, 10 clips/cartridge, 18 cartridges/box

**WK533735 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small w/tape, 25 clips/cartridge, 12 cartridges/box

**WK533835 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small w/tape, 10 clips/cartridge, 18 cartridges/box

**WK534737 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small Strongpoint, 25 clips/cartridge, 12 cartridges/box

**WK534837 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small Strongpoint, 10 clips/cartridge, 18 cartridges/box

**WK534735 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS**

Small Strongpoint w/tape, 25 clips/cartridge, 12 cartridges/box



WK534835 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Small Strongpoint w/tape, 10 clips/cartridge, 18 cartridges/box



WK533702 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium, 25 clips/cartridge, 10 cartridges/box



WK533802 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium, 10 clips/cartridge, 18 cartridges/box



WK533700 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium w/tape, 25 clips/cartridge, 10 cartridges/box



WK533800 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium w/tape, 10 clips/cartridge, 18 cartridges/box



WK533862 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium Large, 10 clips/cartridge, 12 cartridges/box



WK533860 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Medium Large w/tape, 10 clips/cartridge, 12 cartridges/box

WK533872 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Large, 10 clips/cartridge, 12 cartridges/box



WK533870 : HEMOCLIP PLUS® TITANIUM LIGATING CLIPS

Large w/tape, 10 clips/cartridge, 12 cartridges/box

