

CATALOG # :**PG351675****PILLING® MULTIPURPOSE RETRACTOR COMPLETE****Description:**

A versatile, self-retaining Henly-type retractor designed for use in the following surgical specialties: peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions. Spinal surgery in the cervical and lumbar regions. Orthopedic surgery in shoulder and hip procedures. General surgery for inguinal hernia, appendectomy and many minor procedures.

Pilling® Multipurpose Retractor, complete with:

- 1 Frame
- 3 Center Blades in sizes 3/4" (1.9 cm) 1" (2.5 cm), and 1-1/4" (3.1 cm) deep x 5/8" (1.6 cm) wide
- 3 pair of Side Blades 1" (2.5 cm), 2" (5.1 cm), and 3" (7.5 cm) deep x 7/8" (2.2 cm) wide
- Self-retaining, with interchangeable center and side blades. Ball-snap attachment of side blades. Swivel side blades, rotate 360 degrees.

