

DILATORS, BOUGIES, SOUNDS, PROBES, GUIDES**PG355091 : GARRETT VASCULAR DILATOR**

Long 1.0 mm, malleable shaft, 8-1/4" (21 cm)

**PG355092 : GARRETT VASCULAR DILATOR**

Long 1.5 mm, malleable shaft, 8-1/4" (21 cm)

**PG355093 : GARRETT VASCULAR DILATOR**

Long 2.0 mm, malleable shaft, 8-1/4" (21 cm)

**PG355094 : GARRETT VASCULAR DILATOR**

Long 2.5 mm, malleable shaft, 8-1/4" (21 cm)

**PG355095 : GARRETT VASCULAR DILATOR**

Long 3.0 mm, malleable shaft, 8-1/4" (21 cm)

**PG355096 : GARRETT VASCULAR DILATOR**

Long 3.5 mm, malleable shaft, 8-1/4" (21 cm)



PG355097 : GARRETT VASCULAR DILATOR

Long 4.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355098 : GARRETT VASCULAR DILATOR

Long 4.5 mm, malleable shaft, 8-1/4" (21 cm)



PG355099 : GARRETT VASCULAR DILATOR

Long 5.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355100 : GARRETT DILATOR HOLDER

8-1/8" (20.5 cm)



PG355101 : GARRETT DILATOR HOLDER

5-1/8" (13 cm)



PG355275 : HIEBERT DILATOR

0.5mm x 1mm, Malleable shafts, 6" (15.5cm)



PG355276 : HIEBERT DILATOR

1.5mm x 2mm, Malleable shafts, 6" (15.5cm)



PG355277 : HIEBERT DILATOR

2.5mm x 3mm, Malleable shafts, 6" (15.5cm)



PG355278 : HIEBERT DILATOR

3.5mm x 4mm, Malleable shafts, 6" (15.5cm)



PG355285 : DEBAKEY DILATORS 12/SET

19cm




PG355287 : DEBAKEY DILATORS SET

Without case, 19cm




PG355288 : DEBAKEY DILATORS 5/SET

1.0mm, 1.5mm, 2.0mm, 2.5mm, 3.0mm, 19cm



PG355290 : CORONARY SET OF 12 DEBAKEY DILATORS

15" (38cm)




PG355292 : CORONARY SET OF 12 DEBAKEY DILATORS (W/O CASE)

15" (38cm)




AVAILABLE




PG355295 : DEBAKEY DILATOR

1.0mm, 7 1/2" (19cm)




PG355296 : DEBAKEY DILATOR

1.5mm, 7 1/2" (19cm)




PG355297 : DEBAKEY DILATOR

2.0mm, 7 1/2" (19cm)



PG355298 : DEBAKEY DILATOR

2.5mm, 7 1/2" (19cm)



PG355299 : DEBAKEY DILATOR

3.0mm, 7 1/2" (19cm)