


FORCEPS, HEMOSTATS, CLAMPS**PG342611 : MACQUIGG FORCEPS**

Curved, Longitudinal Serrations, 1 5/8" (4cm) tip length, 8 5/8" (22cm)

**PG342615 : FINOCHIETTO THORACIC FORCEPS**

Angled, Longitudinal Serrations, 1 1/8" (2.7cm) tip length, 9 1/2" (24cm)

**PG342625 : RIENHOFF SWAN NECK CLAMP**


Serrated, Curved shaft, 10 5/8" (27cm)

**PG342632 : HARRINGTON CLAMP**


Heavy construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)

**PG342633 : HARRINGTON CLAMP**

Light construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)


**PG342660 : COLLINS LUNG TISSUE FORCEPS 1/2" WIDE**

DeBakey type jaw, 8" (21 cm)




PG342662 : COLLINS LUNG TISSUE FORCEPS 1" WIDE

DeBakey type jaw, 8" (21 cm)




PG342664 : COLLINS LUNG TISSUE FORCEPS 1" WIDE

DeBakey type jaw, 9" (23 cm)




PG342670 : DUVAL LUNG TISSUE FORCEPS 1" WIDE

DeBakey type jaw, 8" (21 cm)



PG342672 : DUVAL LUNG TISSUE FORCEPS 1" WIDE

DeBakey type jaw, 9" (23 cm)



PG342675 : LOVELACE LUNG TISSUE FORCEPS 1" WIDE

Angled on Side, DeBakey type jaw, 7 1/2" (19 cm)




PG342710 : SEMB DISSECTING FORCEPS

Lightly Curved, Smooth, 10 1/4" (26 cm)




PG342711 : SEMB LIGATURE FORCEPS

Acutely Curved, DeBakey type jaw, 9 3/4" (25 cm)



PG342720 : WILLAUER ALLIS TISSUE FORCEPS

Straight, 5 x 6 Teeth, 10" (26 cm)




PG342721 : WILLAUER ALLIS TISSUE FORCEPS

Straight, 5 x 6 Teeth, 11 1/4" (29 cm)




PG342725 : BABCOCK TISSUE FORCEPS

Serrated, 9 1/2" (24 cm)




PG342726 : BABCOCK TISSUE FORCEPS

Serrated, 9 1/2" (24 cm)




PG342800 : LEES BRONCHUS CLAMP

Right Angle, DeBakey Jaw, 1 3/4" (4.5cm) tip length, 8 1/8" (20.5cm)



PG342802 : LEES BRONCHUS CLAMP

Angled, DeBakey Jaw, 1 5/8" (4cm) tip length, Curved shaft, 10" (25.5cm)



PG342804 : BRONCHUS CLAMP

Angled, DeBakey Jaw, 3 1/4" (8.3cm) tip length, 9" (23cm)



PG342822 : SAROT BRONCHUS CLAMP

8.9cm tip length, 9 1/4" (23.5cm)



PG342823 : SAROT BRONCHUS CLAMP

3 1/2" 8.9cm tip length, 9 1/4" (23.5cm)



PG342910 : C.W.MAYO TISSUE FORCEPS

Straight, Serrated, 6.7mm tip, 9" (23 cm)



PG342915 : DAVIS TISSUE FORCEPS

Straight, Serrated, 7.5mm tip, 10" (26 cm)



PG342933 : SINGLEY TISSUE FORCEPS

Serrated, 8.5mm ring end, 9" (23 cm)

