

**FORCEPS, HEMOSTATS, CLAMPS****PG342935 : TUTTLE TISSUE FORCEPS**

Serrated, 6.8mm ring tips, 9" (23 cm)

**PG342950 : NELSON TISSUE FORCEPS**

6 x 7 Teeth, 9" (23 cm)

**PG342965 : BROWN TISSUE FORCEPS**

9 x 9 Teeth, 8" (21 cm)

**PG342967 : BROWN TISSUE FORCEPS**

9 x 9 Teeth, 9 1/2" (24 cm)

**PG342980 : DEBAKEY THORACIC TISSUE FORCEPS**

Straight, 3mm tip, 6 1/8" (16 cm)

**PG342981 : DEBAKEY THORACIC TISSUE FORCEPS**

Straight, 3mm tip, 7 3/4" (20 cm)



**PG342982 : DEBAKEY THORACIC TISSUE FORCEPS**

Straight, 3mm tip, 9 1/2" (25 cm)



**PG342984 : DEBAKEY THORACIC TISSUE FORCEPS**

Straight, 3mm tip, 12" (31 cm)




**PG342986 : DEBAKEY THORACIC TISSUE FORCEPS**

Angled, 3mm tip, 12" (31 cm)



**PG343590 : CHEST DRAIN PULLING CLAMP**

Straight, Serrated, 2" (5cm) tip length, 2mm tip width, 8 5/8" (22cm)



**PG351351 : MORSE TUBING HOLDER FOR SIZE 3/16" AND 1/4" TUBING**

5 3/8" (13.5cm)



**PG351352 : MORSE TUBING HOLDER FOR 5/16" AND 3/8" TUBING**

5 3/8" (13.5cm)



**PG351353 : MORSE TUBING HOLDER FOR 7/16" AND 1/2" TUBING**

5 3/8" (13.5cm)



**PG351354 : MORSE TUBING HOLDER FOR 9/16" AND 5/8" TUBING**

5 3/8" (13.5cm)




**PG351355 : MORSE TUBING HOLDER FOR 11/16" AND 3/4" TUBING**

5 3/8" (13.5cm)



**PG351358 : PRESBYTERIAN HOSPITAL TUBING CLAMP**

Straight, 3/8" (1cm) tip length, 6 1/2" (16.5cm)



**PG351370 : NICHOLS TUBING HOLDER**

7 1/4" (18.5cm)



**PG351373 : VORSE TUBING CLAMP**

Straight, 1 3/8" (3.6cm) tip length, 6 1/8" (15.5cm)



**PG351374 : VORSE TUBING CLAMP**

Straight, 1 5/8" (4cm) tip length, 7 3/8" (19cm)



**PG351375 : VORSE TUBING CLAMP**

Straight, 2" (5.1cm) tip length, 8" (20.5cm)



**PG351376 : VORSE TUBING CLAMP W/GUARD**

Straight, 6 1/8" (15.5cm)



**PG351377 : VORSE TUBING CLAMP W/GUARD**

Straight, 7 3/8" (19cm)



**PG351378 : VORSE TUBING CLAMP W/GUARD**

Straight, 8 1/8" (20.5cm)



**PG351390 : GUNNAR ROLLER CLAMP**

Straight, 7 1/8" (18cm)



**PG351650 : COOLEY VENA CAVA CLAMP**

Curved, 28 Fr. ring tip width, 7 1/4" (18.5cm)

