



FORCEPS, HEMOSTATS, CLAMPS**PG351651 : COOLEY VENA CAVA CLAMP**


Curved, 30 Fr. ring tip width, 7 1/4" (18.5cm)

**PG351653 : COOLEY VENA CAVA CLAMP**


Curved, 36 Fr. ring tip width, 7 1/4" (18.5cm)

**PG351654 : COOLEY VENA CAVA CLAMP**

Curved, 34 Fr. ring tip width, 7 1/4" (18.5cm)

**PG351743 : REVUELTA-GARCIA ALLIS CLAMP**

Half Curved, DeBakey Jaw, 1" (2.4cm) tip length, 10 1/8" (25.5cm)

**PG351744 : REVUELTA-GARCIA ALLIS CLAMP**

Half Curved, DeBakey Jaw, 10 7/8" (27.5cm)

**PG351748 : CARDIAC VALVE FORCEPS**

Serrated, 3mm tip, 9 1/2" (25 cm)




PG351800 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 6" (16 cm)



PG351801 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 6" (16 cm)



PG351802 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 7 3/4" (20 cm)




PG351803 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 7 3/4" (20 cm)



PG351803DG : DEBAKEY-DIETHRICH DIAMA-GRIP TISSUE FORCEPS

Straight, 1.3 mm tip, 7-3/4" (20 cm), Diama-Grip jaws



PG351804 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 9 1/2" (24 cm)



PG351806 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 9 1/2" (24 cm)

PG351808 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 12" (31 cm)



PG351810 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 6" (16 cm)



PG351812 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 7 3/4" (20 cm)



PG351814 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 9 1/2" (24 cm)



PG351815 : DEBAKEY TISSUE FORCEPS

2mm tip, 6" (15 cm), TC



PG351817 : DEBAKEY TISSUE FORCEPS

2mm tip, 7 3/4" (20 cm), TC



PG351818 : DEBAKEY TISSUE FORCEPS

Delicate, 7 3/4" (20 cm), TC



PG351819 : DEBAKEY TISSUE FORCEPS

9 1/2" (24 cm), TC




PG351820 : DEBAKEY ULTRALIGHT TISSUE FORCEPS

1.5mm tip, 6" (16 cm)



PG351821 : DEBAKEY ULTRALIGHT TISSUE FORCEPS

1.5mm tip, 7 3/4" (20 cm)



PG351824 : DEBAKEY TISSUE FORCEPS

Delicate, 9 1/2" (24 cm), TC



PG351825 : WHEELER VESSEL FORCEPS

DeBakey type jaw, 1.6mm , 5 1/2" (14 cm)

