

**FORCEPS, HEMOSTATS, CLAMPS****PG351898 : CRAFOORD TISSUE FORCEPS**

Serrated, 8" (21 cm)

**PG351900 : GERALD DRESSING FORCEPS**

Serrated, 7" (18 cm)

**PG351901 : GERALD ATRAUGRIP TISSUE FORCEPS**

DeBakey type jaw, 7" (18 cm)

**PG351901DG : GERALD DIAMA-GRIP DRESSING FORCEPS**

Straight, 1.1 mm tip, 7" (18 cm), Diama-Grip jaws

**PG351902 : GERALD DRESSING FORCEPS**

Serrated, 9 3/4" (25 cm)

**PG351905 : GERALD TISSUE FORCEPS**

1 x 2 Teeth, 7" (18 cm)



**PG351907 : GERALD TISSUE FORCEPS**

1 x 2 Teeth, 9 3/4" (25 cm)



**PG351908 : GERALD ATRAUGRIP TISSUE FORCEPS**

DeBakey type jaw, 8" (21 cm)



**PG351909 : GERALD ATRAUGRIP TISSUE FORCEPS**

DeBakey type jaw, 9 1/2" (24 cm)



**PG351910 : GERALD ATRAUGRIP TISSUE FORCEPS**

Angled, DeBakey type jaw, 6 7/8" (18 cm)



**PG353001 : JOHNS HOPKINS BULLDOG CLAMP**

Straight, Serrated, 1 1/2" (4cm)



**PG353002 : JOHNS HOPKINS BULLDOG CLAMP**

Straight, Serrated, 1 7/8" (5cm)



**PG353003 : JOHNS HOPKINS BULLDOG CLAMP**

Straight, Serrated, 2 1/4" (5.5cm)

**PG353011 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 1 1/2" (4cm)



**PG353012 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 1 7/8" (5cm)



**PG353014 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 2 1/2" (6.5cm)



**PG353016 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 3" (7.5cm)



**PG353025 : HESS BULLDOG CLAMP**

Angled, Serrated, 1 1/8" (2.7cm) tip length, 2 1/8" (5.5cm)



**PG353035 : BULLDOG CLAMP APPLYING FORCEPS**

Angled, Smooth, 9 1/2" (24 cm)



**PG353040 : SCHWARTZ CLAMP HOLDING FORCEPS**

Angled, 7 1/2" (19 cm)





**PG353041 : SCHWARTZ MICRO BULLDOG CLAMP**

Straight, 1" (2.5cm)



**PG353042 : SCHWARTZ MICRO BULLDOG CLAMP**

Slight Curved, 1/4" (6.4mm) tip length, 1" (2.5cm)



**PG353043 : SCHWARTZ MICRO BULLDOG CLAMP**

Strongly Curved, 3/8" (10.2mm) tip length, 1" (2.5cm)



**PG353050 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 3/4" (2cm) tip length, 3 1/8" (8cm)



**PG353052 : DEBAKEY CROSS ACTION BULLDOG CLAMP**

Straight, 1 1/8" (3cm) tip length, 3 1/2" (9cm)