


**SUCTION, ASPIRATING TUBES, NEEDLES****PG351577 : DEBAKEY HEPARIN NEEDLE LARGE**

Malleable, 4 mm tip, 6" (15.1 cm)

**PG351578 : IMA HEPARIN NEEDLE**

23 gauge tip, 2-1/8" (5.3 cm)

**PG351586 : GARRETT HEPARIN INJECTOR**

Malleable, 5/16" (8 mm) tip, working length 3" (7.5 cm)

**PG351587 : GARRETT HEPARIN INJECTOR**

Malleable, 1/2" (12 mm) tip, working length 3" (7.5 cm)

**PG351593 : SAMSON HEPARIN INJECTOR**

Malleable, small 7 mm tip, 3-3/8" (8.5 cm)

**PG351594 : SAMSON HEPARIN INJECTOR MEDIUM**

Malleable, medium 9 mm tip, 3-3/8" (8.5 cm)



**PG507453 : ASPIRATING NEEDLE**

Reinforced angular shaft, with protective cover, distal needle 23 gauge (0.6 mm), 14.7 mm long, working length 8 1/4" (21 cm)



**PG507454 : SPARE DISTAL NEEDLE TIP FOR ASPIRATING NEEDLE 507453**

23 gauge (0.6 mm), 14.7 mm long



**PG536220 : ANDREWS PYNCHON SUCTION TUBE**

Curved, 6mm tip, 3mm dia. removable shaft, 9 5/8" (24.5cm)



**PO271001 : SUCTION TUBES SANTY L=28CM**




**PO271002 : SUCTION TUBES SANTY L=26,5CM**



**PO271017 : HEPARIN CANNULA N°1**



**PO271018 : HEPARIN CANNULA N° 2**



**PO271019 : HEPARIN CANNULA N° 3**



**PO271020 : HEPARIN CANNULA FINE MODEL**



**PO271021 : SUCTION TUBE NOWACK 23CM DIAM 4MM**



**PO271052 : SUCTION TUBES CONE-LUER ANGLED D=0,9MM**



**PO271053 : SUCTION TUBE CONE LUER 10CM ANGLED DIAM 1.2MM**



**PO271055 : SUCTION TUBES CONE-LUER ANGLED L=11,5CM D=1,8MM**



**PO271056 : SUCTION TUBES CONE-LUER ANGLED L=10CM D=1,9MM**



**PO271057 : SUCTION TUBE CONE LUER 10CM ANGLED DIAM  
2.4MM**



**PO271067 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM  
D=0,9MM**



**PO271068 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM  
D=1,5MM**



**PO271069 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM  
D=1,9MM**



**PO271070 : SUCTION TUBE 12CM CURVED OLIVE SHAPED TIP  
2.5MM**

