


NEEDLE HOLDERS

**PG343000 : GERBODE NEEDLE HOLDER**

16.5 x 2mm jaw, 10 1/2" (26.5cm), TC

**PG343027 : STRATTE NEEDLE HOLDER**

Reverse curved, 9" (23cm), TC

**PG343030 : BERRY STERNAL "BIG UGLY" TECHNIQUE® NEEDLE HOLDER**

heavy-duty needle holder with 6 mm-wide jaws designed to securely hold the large needles used to close sternotomies. LIFE-LOK Box Lock.

**PG354900 : DEBAKEY NEEDLE HOLDER**


13.7 x 2 mm jaw, 7" (18 cm), TC

**PG354902 : DEBAKEY NEEDLE HOLDER**

17 x 1.8 mm jaw, 9-1/4" (23.5 cm), TC


**PG354904 : DEBAKEY NEEDLE HOLDER**

21.6 x 2.3mm jaw, 10 1/2" (26.5cm), TC




PG354906 : DEBAKEY NEEDLE HOLDER

19 x 2.6mm jaw, 11 3/4" (30cm), TC




PG354915 : W.M.LEMMON NEEDLE HOLDER

1.5mm jaw, 5 3/4" (14.5cm), TC



PG354916 : W.M.LEMMON NEEDLE HOLDER

1.5mm jaw, 6 3/4" (17cm), TC




PG354917 : W.M. LEMMON NEEDLE HOLDER

7.2 x 1.6mm jaw, 7 1/2" (19cm), TC




PG354920 : HOWELL CORONARY NEEDLE HOLDER

14 x 1.6mm jaw, 7 1/4" (18.5cm), TC




PG354928 : HOWELL VALVE NEEDLE HOLDER

15.3 x 2mm jaw, 8" (20.5cm), TC



PG354931 : BAILEY NEEDLE HOLDER

17 x 2mm jaw, 7" (18.5cm), TC



PG354941 : JAMISON NEEDLE HOLDER

1.8mm jaw, 9" (23cm), TC



PG354958 : PILLING THIN VASCULAR NEEDLE HOLDER

16.3 x 1.8mm jaw, 9" (23cm), TC



PG354959 : PILLING THIN VASCULAR NEEDLE HOLDER

16 x 2.2mm jaw, 10 1/4" (26cm), TC



PG354960 : PILLING THIN VASCULAR NEEDLE HOLDER

11 x 1.8mm jaw, 6" (15cm), TC



PG354961 : PILLING THIN VASCULAR NEEDLE HOLDER

11.8 x 1.7mm jaw, 7" (17.5cm), TC



PG354962 : PILLING THIN VASCULAR NEEDLE HOLDER

15.8 x 1.8 mm jaw, 8-1/4" (21 cm), TC



PG354964 : CASTROVIEJO TECHNIQUE® NEEDLE HOLDERS

Straight, with lock, 13 x 1.4 mm serrated jaw, 7-1/4" (18.5 cm), TC



PG354964DG : CASTROVIEJO NEEDLE HOLDER

Straight, Flat handle, With lock, 1.8mm Diama-Grip jaw, 7 1/4" (18.5cm)



PG354965 : CASTROVIEJO TECHNIQUE® NEEDLE HOLDERS

Straight, with lock, 13 x 1.4 mm serrated jaw, 8" (20.5 cm), TC



PG354965DG : CASTROVIEJO NEEDLE HOLDER

Flat handle, w/Lock, Diama-Grip, 8" (20.5cm)



PG354966 : MICRO SURGERY TECHNIQUE® NEEDLE HOLDERS

Straight, Flat handle, With lock, 13 x 1.4mm smooth jaw, 7" (18cm), TC



PG354967 : MICRO SURGERY TECHNIQUE® NEEDLE HOLDERS

Straight, flat handle, with lock, 13 x 1.4 mm smooth jaw, 8" (20.5 cm), TC