

RETRACTORS, HOOKS**PO161073 : RET OBWEGESER FOR SOFT TISSUES CURVED DOWN
22CM / 80MMX16MM****PO161074 : RET OBWEGESER FOR SOFT TISSUES CURVED DOWN
22CM / 70MMX14MM****PO161075 : RET OBWEGESER FOR SOFT TISSUES CURVED DOWN
22CM / 55MMX12MM****PO161076 : RET OBWEGESER FOR SOFT TISSUES CURVED DOWN
22CM / 42MMX11MM****PO161077 : RET OBWEGESER FOR SOFT TISSUES CURVED DOWN
22CM / 31MMX10MM****PO161078 : RETRACTOR OBWEGESER FOR SOFT TISSUES
CURVED DOWN 25X7MM**



**PO161080 : RET OBWEGESER FOR SOFT TISSUES CUTTING TOP
22CM / 80MMX16MM**



**PO161081 : RET OBWEGESER FOR SOFT TISSUES CUTTING TOP
22CM / 60MMX12MM**



**PO161082 : RET OBWEGESER FOR SOFT TISSUES CUTTING TOP
22CM / 42MMX13MM**




**PO161092 : RETRACTORS WITH HANDLE MICKULITZ L=27,5CM
W=85MM H=35MM**




PO161104 : RETRACTOR WITH HANDLE KOCHER 22CM/35X8MM



PO161105 : RETRACTOR WITH HANDLE KOCHER 22CM/40X11MM



PO161106 : RETRACTOR WITH HANDLE KOCHER 22CM/55X11MM



PO161107 : RETRACTOR WITH HANDLE KOCHER 22CM/70X14MM



PO161108 : RETRACTOR DAUTREY SMALL MODEL 14CM



PO161109 : RETRACTOR WITH HANDLE DAUTREY LARGE MODEL 14CM



PO161112 : RETRACTOR CORTY 20CM



PO164004 : RETRACTOR SELF RETAINING BACKMANN 13CM FOR HAND



PO164005 : RETRACTORS SELF-RETAINING BECKMANN L=14CM



PO164006 : RETRACTORS SELF-RETAINING BECKMANN L=17CM



PO164007 : RETRACTORS SELF-RETAINING BECKMANN L=20,5CM



PO164008 : RETRACTORS SELF-RETAINING BECKMANN PRONGS L=31CM W=20MM



PO164009 : RETRACTORS SELF-RETAINING BECKMANN PRONGS L=31CM W=30MM



PO164024 : RETRACTOR BECKMANN STILLE 33CM DEPTH 35MM BLUNT PRONGS



PO164025 : RETR SELF-RETAIN BECKMANN-STILLE SHARP PRONGS L=33CM W=35MM

