

**RETRACTORS, HOOKS****PG164543 : RICHARDSON RETRACTOR**

1" (2.5cm) width x 1 1/4" depth blade, 9 3/4" (25.0cm)

**PG164544 : RICHARDSON RETRACTOR**

1 1/2" (3.8cm) width x 1 1/2" depth blade, 9 3/4" (25.0cm)

**PG164550 : RICHARDSON-EASTMAN RETRACTOR**

Set of 2. Large 10 3/4", 65mm x 43mm and 42mm x 38mm blades, Small 10" (25.5cm), 33mm x 30mm and 26mm x 23mm blades.

**PG164630 : LANGENBECK RETRACTOR**

3/8" (1.0cm) width x 1 5/8" (4.0cm) depth blade, 8 1/2" (21.5cm)

**PG164631 : LANGENBECK RETRACTOR**

1/2" (1.3cm) width x 1 5/8" (4.0cm) depth blade, 8 1/2" (21.5cm)

**PG164632 : LANGENBECK RETRACTOR**

5/8" (1.6cm) width x 1 5/8" (4.0cm) depth blade, 8 1/2" (21.5cm)



**PG164635 : SAUERBRUCH RETRACTOR**

5/8" (1.5cm) width x 1 3/4" (4.5cm) depth blade, 9" (23.0cm)



**PG164636 : SAUERBRUCH RETRACTOR**

7/8" (2.3cm) width x 2 1/2" (6.3cm) depth blade, 9" (23.0cm)



**PG164637 : SAUERBRUCH RETRACTOR**

3/4" (2.0cm) width x 2 7/8" (7.3cm) depth blade, 8 7/8" (22.5cm)



**PG164650 : BREWSTER RETRACTOR**

3/4" (2.0cm) width x 1 3/4" (4.5cm) depth blade, 8 1/2" (21.5cm)




**PG164652 : BREWSTER RETRACTOR**

1" (2.5cm) width x 2" (5.0cm) depth blade, 9" (23.0cm)



**PG164656 : BREWSTER RETRACTOR**

1 3/4" (4.5cm) width x 2 1/2" (6.5cm) depth blade, 9 1/4" (23.5cm)




**PG164658 : BREWSTER RETRACTOR**

1 7/8" (4.7cm) width x 2 7/8" (7.0cm) depth blade, 9 1/2" (24.0cm)



**PG164698 : PARKER-MOTT RETRACTOR**

1" (2.4cm) / 1" (2.4cm) width x 1 1/4" (3.0cm) / 2 1/2" (65.0cm) depth blade, 5 3/4" (14.5cm), 2/set




**PG164700 : FARABEUFR RETRACTOR**

5/8" (1.5cm) / 5/8" (1.5cm) width x 1 1/8" (2.5cm) / 1 1/4" (3.0cm) depth blade, 6" (15.0cm), 2/set



**PG164702 : FARABEUFR RETRACTOR**

3/8" (1.1cm) / 1/2" (1.2cm) width blade, 4 1/2" (11.5cm), 2/set



**PG164710 : MAYO-COLLINS RETRACTOR**

3/4" (1.9cm) / 3/4" (1.9cm) width x 1 1/8" (3.0cm) / 1 3/8" (3.5cm) depth blade, 6" (15.0cm), 2/set



**PG164715 : ARMY-NAVY RETRACTOR**

5/8" (1.5cm) width blade, 8 1/2" (21.5cm), 2/set




**PG164750 : SENN RETRACTOR**

Sharp, 3-pronged, 1/4" (6mm) width x 3/4" (2.0cm) depth blade, 6 1/4" (16.0cm)



**PG164752 : SENN RETRACTOR**

Blunt, 3-pronged, 1/4" (6mm) width x 3/4" (2.0cm) depth blade, 6 1/4" (16.0cm)




**PG164754 : CHILDREN'S RETRACTOR**

Blades 20mm long x 9.5mm wide (3/4 x 3/8) / 20mm long x 6mm wide (3/4 x 1/4), 6 1/8" (15.5cm)



**PG164790 : MEYERDING FINGER RETRACTOR**

1/4" (7.5mm) width blade, 6 7/8" (17.5cm)



**PG164865 : CUSHING VEIN RETRACTOR**

1/2" (11.8mm) width x 1/2" (1.0cm) depth blade, 8 1/2" (21.5cm)



**PG165085 : HEISS RETRACTOR**

4 x 4 Sharp prongs, Side Blade 3/8" (9.5mm) wide, 4 1/8" (10.5cm)



**PG165086 : HEISS RETRACTOR**

4 x 4 Blunt prongs, Side Blade 3/8" (9.5mm) wide, 4 1/8" (10.5cm)