

RETRACTORS, HOOKS



PG341390 : TUFFIER RETRACTOR

Spread 5 7/8" (15cm), Scored blades, Side Blade 1 13/16" (4.6cm) wide x 2 3/16" (5.6cm) deep, Arm length 4 1/2" (11.4cm)



PG341397 : KUYPER-MURPHY STERNAL RETRACTOR, SMALL

Concave blade, Spread 4 1/2" (11.5cm), Side Blade 3 1/2" (9cm) wide x 1 1/8" (3cm) deep, Straight arms



PG341398 : KUYPER-MURPHY STERNAL RETRACTOR, MEDIUM

Concave blade, Spread 6 5/16" (16cm), Side Blade 3 1/2" (9cm) wide x 1 1/8" (3cm) deep, Straight arms



PG341448 : BAILEY RIB APPROXIMATOR RETRACTOR, SMALL

Sharp




PG341450 : BAILEY RIB APPROXIMATOR RETRACTOR, ADULT

Sharp




PG341460 : BAILEY-GIBBON RIB APPROXIMATOR RETRACTOR

Sharp




PG341501 : LEMOLE-PILLING VALVE RETRACTOR

Set #1 fits Cooley and Chaux Sternal Retractors



PG341502 : LEMOLE-PILLING VALVE RETRACTOR

Set #2 fits Morse-style Retractors



PG341503 : LEMOLE-PILLING VALVE RETRACTOR

Set #3 fits Ankeney Retractor



PG341506 : LEMOLE VALVE RETRACTOR #1 INCLUDES 341195 COOLEY RETRACTOR



PG341511 : LEMOLE-PILLING T FRAME ONLY FOR TYPE #1 RETRACTOR




PG341512 : LEMOLE-PILLING T FRAME ONLY FOR TYPE #2 RETRACTOR



PG341513 : LEMOLE-PILLING T FRAME ONLY FOR TYPE #3 RETRACTOR

PG341514 : CHAUX-LEMOLE BI-ANGLE T-FRAME FOR COOLEY AND CHAUX RETRACTORS



PG341515 : LEMOLE SIDE ARM EXTENSION FOR COOLEY AND CHAUX RETRACTORS



PG341520 : UNIVIVERSAL CLAMP FOR LEMOLE-PILLING RETRACTORS




PG341522 : COOLEY ATRIAL RETRACTOR, SMALL

5/8" (1.5cm) width blade, 9 1/2" (24.0cm)




PG341523 : COOLEY ATRIAL RETRACTOR, MEDIUM

7/8" (2.0cm) width blade, 9 1/2" (24.0cm)



PG341524 : COOLEY ATRIAL RETRACTOR, LARGE



1 3/8" (3.5cm) width blade, 9 1/2" (24.0cm)



PG341525 : MALLEABLE ATRIAL RETRACTOR

3/4" (2.0cm) width blade, 8 3/4" (22.0cm)





PG341635 : BREWSTER RETRACTOR

10mm / 13mm & 10mm / 13mm width blade



PG341640 : SEMB RETRACTOR

3/4" (2.0cm) width x 1 3/4" (4.5cm) depth blade, 9 3/4" (25.0cm)



PG341650 : DAVIDSON SCAPULA RETRACTOR, SMALL

2" (5cm) width blade, 7 3/4" (19.5cm)



PG341651 : DAVIDSON SCAPULA RETRACTOR, LARGE

3 1/2" (9cm) width blade, 7 3/4" (19.5cm)



PG341678 : ALLISON LUNG RETRACTOR, CHILD

1 1/2" (4.0cm) width blade, 10 1/4" (26.0cm)