

RETRACTORS, HOOKS**PG351669 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR**

Sharp, Center blade 5/8" (1.6cm) wide x 7/8" (2.2cm) deep

**PG351670 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR**

Sharp, Center blade 5/8" (1.6cm) wide x 1 1/8" (2.8cm) deep

**PG351671 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR**

Sharp, Center blade 5/8" (1.6cm) wide x 1 1/4" (3.2cm) deep

**PG351672 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR**

Sharp, Center blade 5/8" (1.6cm) wide x 2" (5.2cm) deep

**PG351673 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR**

Sharp, Center blade 5/8" (1.6cm) wide x 2 1/2" (6.5cm) deep

**PG351675 : PILLING® MULTIPURPOSE RETRACTOR COMPLETE**

A versatile, self-retaining Henly-type retractor designed for use in the following surgical specialties: peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions. Spinal surgery in the...



PG351676 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 1" (2.5cm) deep, 1 pair



PG351677 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 2" (5cm) deep, 1 pair



PG351678 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 3" (7.6cm) deep, 1 pair




PG351679 : IMPROVED HENLY RETRACTOR WITH STRAIGHT ARMS

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



PG351680 : IMPROVED HENLY RETRACTOR WITH ARMS SPREAD OUTWARD

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



PG351747 : REVUELTA GARCIA VALVE LEAFLET RETRACTOR

Curved, 24.0mm width blade



PG351750 : BLANCO VALVE RETRACTOR

5mm width x 3/8" (1.0cm) depth blade, 11 1/2" (29.0cm)

PG351752 : BLANCO VALVE RETRACTOR

7mm width x 3/8" (1.0cm) depth blade, 11" (28.0cm)




PG351754 : CUSHING-BLANCO VALVE RETRACTOR

12mm width x 1/4" (0.5cm) depth blade, 11" (28.0cm)



PG351755 : BLANCO VALVE RETRACTOR

9mm width x 1/2" (1.0cm) depth blade, 11 1/2" (29.0cm)



PG351763 : ABLAZA-BLANCO RETRACTOR

9/16" (1.4cm) width blade, 13" (33.0cm)



PG351767 : DESMARRE VALVE RETRACTOR

4mm width blade, 7 1/4" (18.5cm)



PG351768 : DESMARRE VALVE RETRACTOR

6mm width blade, 7 1/4" (18.5cm)



PG351769 : DESMARRE VALVE RETRACTOR

8mm width blade, 7 1/4" (18.5cm)



PG351770 : DESMARRE VALVE RETRACTOR

10mm width blade, 7 1/4" (18.5cm)



PG351771 : DESMARRE VALVE RETRACTOR

12mm width blade, 7 1/4" (18.5cm)



PG351772 : DESMARRE VALVE RETRACTOR

15mm width blade, 7 1/4" (18.5cm)



PG351773 : CRAWFORD AORTIC RETRACTOR

15mm width x 3/4" (2.0cm) depth blade, 7 3/8" (18.5cm)



PG351774 : CRAWFORD AORTIC RETRACTOR

22mm width x 3/4" (2.0cm) depth blade, 7 7/8" (20.0cm)

