

CVT



PG359187 : DEBAKEY TITANIUM OCCLUSION CLAMP

Angled, 9 5/8" (24.5cm)




PG359188 : DEBAKEY TITANIUM OCCLUSION CLAMP

Angled, 1 3/4" (4.6cm) tip length, 9 1/2" (24cm)




PG359191 : SATINSKY TITANIUM CLAMP

Angled, DeBakey Jaw, 1 3/4" (4.45cm) tip length, 9 5/8" (24.5cm)



PG359193 : SATINSKY TITANIUM CLAMP

Angled, DeBakey Jaw, 9 1/2" (24cm)




PG359194 : GLOVER TITANIUM ARICULAR CLAMP

Angled, DeBakey Jaw, 1 5/8" (4.2cm) tip length, 9 3/4" (24.5cm)



PG359195 : DEBAKEY TITANIUM VASCULAR CLAMP

Angled 53, DeBakey Jaw, 10 1/2" (26.5cm)




PG359201 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 10 GRAMS

Straight, DeBakey type jaw, 5 1/8" (13 cm)



PG359206 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 10 GRAMS

Curved, DeBakey type jaw, 4 3/4" (12 cm)



PG359208 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 10 GRAMS

Curved, DeBakey type jaw, 6 1/2" (17 cm)




PG359209 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 10 GRAMS

Curved, DeBakey type jaw, 6 1/2" (17 cm)




PG359210 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 10 GRAMS

Curved, DeBakey type jaw, 6 3/8" (17 cm)



PG359211 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 15 GRAMS

Straight, DeBakey type jaw, 6 5/8" (17 cm)



PG359212 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 15 GRAMS

Angled, DeBakey type jaw, 6 1/4" (16 cm)

PG359213 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 15 GRAMS

Angled, DeBakey type jaw, 6" (15 cm)




PG359214 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 15 GRAMS

Angled, DeBakey type jaw, 5 7/8" (15 cm)




PG359215 : CASTANEDA TITANIUM VASCULAR CLAMPING FORCEPS 15 GRAMS

Angled, DeBakey type jaw, 5 1/2" (14 cm)



PG359221 : COOLEY TITANIUM PEDIATRIC CLAMP

Spoon, Serrated, 2.2mm tip width, 7 1/4" (18.5cm)



PG359236 : DERRA TITANIUM ANASTAMOSIS CLAMP

Curved, DeBakey Jaw, 1 1/8" (2.7cm) tip length, 2.6mm tip width, 6 3/4" (17cm)



PG359237 : DERRA TITANIUM ANASTAMOSIS CLAMP

Curved, DeBakey Jaw, 2.8mm tip width, 6 3/4" (17.5cm)



PG359238 : RENAL ARTERY TITANIUM CLAMP

Angled, Serrated, 2.5mm tip width, 7 1/2" (19cm)





PG359250 : MICRO TITANIUM SCISSORS

Curved, Serrated, 8 1/4" (21.0cm)



PG359253 : MICRO TITANIUM SCISSORS

Straight, 9" (23.0cm)



PG359259 : MICRO TITANIUM SCISSORS

Round Handle, Straight, 6 1/4" (16.0cm)



PG359260 : MICRO TITANIUM SCISSORS

Round Handle, Curved, 6 1/4" (16.0cm)



PG359261 : MICRO TITANIUM SCISSORS

Round Handle, Straight, Serrated, 6 1/4" (16.0cm)

