
**PG841050 : KAISER-PILLING ATRAUGRIP RING FORCEPS (1AG)**



Straight, DeBakey type jaw, 12-5/8" (32 cm)


**PG841051 : KAISER-PILLING ATRAUGRIP RING FORCEPS (5AG)**

Reverse curved, DeBakey type jaw, 11-3/4" (30 cm)


**PG841052 : KAISER-PILLING ATRAUGRIP RING FORCEPS (2AG)**

Curved, DeBakey type jaw, 12-1/4" (31 cm)


**PO354002 : REPLACEMENT STILETTE ADULT****PO354003 : REPLACEMENT STILETTE CHILD****PO354004 : DILATOR VOLLMAN 50CM RING 4MM**




PO354005 : DILATOR VOLLMAN 50CM RING 6MM




PO354006 : DILATOR VOLLMAR 50CM RING 7MM




PO354007 : DILATOR VOLLMAR 50CM RING 8MM



PO354008 : DILATOR VOLLMAR 50CM RING 10MM



PO354009 : DILATOR VOLLMAR 50CM RING 12MM



PO354010 : DILATOR VOLLMAR 50CM RING 14MM



PO354045 : DILATOR DE BAKEY 24CM 1MM

PO354046 : DILATOR DE BAKEY 24CM 2MM



PO354047 : DILATOR DE BAKEY 24CM 3MM



PO354048 : DILATOR DE BAKEY 24CM 4MM



PO354049 : DILATOR DE BAKEY 24CM 5MM




PO354050 : DILATOR DE BAKEY 24CM 6MM



PO354051 : DILATOR DE BAKEY 24CM 7MM



PO354052 : DILATOR DE BAKEY 24CM 8MM





PO354053 : DILATOR DE BAKEY 24CM 9MM



PO354054 : DILATOR DE BAKEY 24CM 10MM



PO355330 : NEEDLE HOLDER CASTROVIEJO CURVED FLAT HANDLE 23CM TC



PO355331 : NEEDLE HOLDER CASTROVIEJO CURVED FLAT HANDLE WITH LOCK 23CM TC



PO355341 : NEEDLE HOLDER CASTROVIEJO STRAIT FLAT HANDLE 23CM TC

