

**PG341696 : HAIGHT RETRACTOR**

4 3/4" (12cm) width blade, 12 3/4" (32.5cm)

**PG341698 : HAIGHT RETRACTOR**

5 1/2" (14cm) width blade, 13 1/2" (34.5cm)

**PG341720 : WILLAUER-DEAVER RETRACTOR**

Hollow handle, 1" (2.5cm) width x 6 1/4" (15.9cm) depth blade, 14 1/8" (36.0cm)

**PG341722 : WILLAUER-DEAVER RETRACTOR**

Hollow handle, 2" (5cm) width x 6 1/4" (15.9cm) depth blade, 14 1/8" (36.0cm)

**PG341730 : LEMMON EXTRA DEEP DEAVER RETRACTOR**

Hollow handle, 1" (2.5cm) width x 9 1/2" (24.3cm) depth blade, 15" (38.0cm)

**PG341732 : LEMMON EXTRA DEEP DEAVER RETRACTOR**

Hollow handle, 2" (5cm) width x 9 1/2" (24.3cm) depth blade, 15" (38.0cm)



PG341740 : DEBAKEY-COOLEY RETRACTOR

Hollow handle, 2 1/4" (5.7cm) width x 7 1/4" (18.4cm) depth blade, 12 1/4" (31.0cm)




PG341770 : HARRINGTON SPLANCHNIC RETRACTOR

Hollow handle, 1 1/2" (3.8cm) width x 4 3/4" (12.0cm) depth blade, 12 1/2" (31.5cm)



PG341773 : HARRINGTON SPLANCHNIC RETRACTOR

Hollow handle, 2 1/2" (6.5cm) width x 5" (12.7cm) depth blade, 12" (30.5cm)



PG341780 : WYLIE SPLANCHNIC RETRACTOR

Hollow handle, 2 1/2" (6.3cm) width x 6 3/4" (17.1cm) depth blade




PG341790 : HARRINGTON EXTRA DEEP RETRACTOR

Hollow handle, 2 1/4" (5.7cm) width x 8" (20.5cm) depth blade, 11 3/4" (30.0cm)



PG341791 : HARRINGTON EXTRA DEEP RETRACTOR

Hollow handle, 4 1/2" (11.5cm) width x 7" (18 cm) depth blade, 13 1/8" (33.5cm)



PG341795 : WYLIE RENAL VEIN RETRACTOR

1" (2.5cm) width x 3 3/4" (9.5cm) depth blade, 11 3/4" (30.0cm)



PG341797 : WYLIE RENAL VEIN RETRACTOR

1" (2.5cm) width x 6 7/8" (17.5 cm) depth blade, 12 3/4" (32.5cm)



PG341803 : PILLING STILLE-LUER DUCKBILL RONGEUR

Double action, Curved on side, 9" (23cm)



PG341836 : PILLING RUSKIN RONGEUR

Double action, Curved, 4mm bite, 6 1/8" (15.5cm)



PG341992 : ROOS BRACHIAL PLEXUS ROOT RETRACTOR

Angled, 13.5mm width blade



PG342002 : WOLVEK STERNAL APPROXIMATOR

10 3/8" (26.5cm)



PG342004 : WOLVEK STERNAL CRIMPER



PG342006 : STERNAL LOCKING PLATES

An implantable device used to fixate size 5 (B&S Gage 20) monofilament stainless steel suture wire (Deknatel 29727, Ethicon M650) to achieve sternal approximation. Only 316L stainless steel sternal closure wire should be...



PG342015 : GRUNWALD SUTURE RING

13 1/4" (33.6cm)




PG342020 : LEBSCHKE STERNUM KNIFE

Straight tip, Straight shaft, 9 7/8" (25cm)



PG342025 : PARSONNET EPICARDAL RETRACTOR

3 X 3-pronged, 1 7/8" (4.5cm)



PG342040 : SWEET MODIFIED STERNAL PUNCH

10" (25.5cm)



PG342043 : KIRKLIN STERNAL AWL

Working length 2" (5.0cm), Pointed needle end, Wrench included, 7 1/8" (18.2cm)

