**PG342601 : MIXER FORCEPS**

Strongly Curved, Serrated, 3/4" (2cm) tip length, 8 3/4" (22cm)

**PG342603 : MIXER FORCEPS**

Half Curved, Serrated, 5/8" (1.65cm) tip length, 2.3mm tip width, 9 1/8" (23cm)

**PG342605 : MIXER FORCEPS**

Strongly Curved, DeBakey Jaw, 1 3/8" (3.5cm) tip length, 11" (28cm)

**PG342606 : MIXER FORCEPS**


Strongly Curved, Serrated (3cm) tip length, 10 3/4" (27.5cm)

**PG342611 : MACQUIGG FORCEPS**

Curved, Longitudinal Serrations, 1 5/8" (4cm) tip length, 8 5/8" (22cm)

**PG342615 : FINOCHIETTO THORACIC FORCEPS**

Angled, Longitudinal Serrations, 1 1/8" (2.7cm) tip length, 9 1/2" (24cm)



**PG342625 : RIENHOFF SWAN NECK CLAMP**

Serrated, Curved shaft, 10 5/8" (27cm)




**PG342632 : HARRINGTON CLAMP**

Heavy construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)




**PG342633 : HARRINGTON CLAMP**

Light construction Right Angle, Longitudinal Serrations, 1 1/2" (3.85cm) tip length, 12 3/8" (31.5cm)




**PG342660 : COLLINS LUNG TISSUE FORCEPS 1/2" WIDE**

DeBakey type jaw, 8" (21 cm)




**PG342662 : COLLINS LUNG TISSUE FORCEPS 1" WIDE**

DeBakey type jaw, 8" (21 cm)




**PG342664 : COLLINS LUNG TISSUE FORCEPS 1" WIDE**

DeBakey type jaw, 9" (23 cm)



**PG342670 : DUVAL LUNG TISSUE FORCEPS 1" WIDE**

DeBakey type jaw, 8" (21 cm)



**PG342672 : DUVAL LUNG TISSUE FORCEPS 1" WIDE**

DeBakey type jaw, 9" (23 cm)




**PG342675 : LOVELACE LUNG TISSUE FORCEPS 1" WIDE**

Angled on Side, DeBakey type jaw, 7 1/2" (19 cm)




**PG342710 : SEMB DISSECTING FORCEPS**

Lightly Curved, Smooth, 10 1/4" (26 cm)



**PG342711 : SEMB LIGATURE FORCEPS**

Acutely Curved, DeBakey type jaw, 9 3/4" (25 cm)




**PG342720 : WILLAUER ALLIS TISSUE FORCEPS**

Straight, 5 x 6 Teeth, 10" (26 cm)



**PG342721 : WILLAUER ALLIS TISSUE FORCEPS**


Straight, 5 x 6 Teeth, 11 1/4" (29 cm)



**PG342725 : BABCOCK TISSUE FORCEPS**


Serrated, 9 1/2" (24 cm)

AVAILABLE




**PG342726 : BABCOCK TISSUE FORCEPS**

Serrated, 9 1/2" (24 cm)




**PG342800 : LEES BRONCHUS CLAMP**

Right Angle, DeBakey Jaw, 1 3/4" (4.5cm) tip length, 8 1/8" (20.5cm)



**PG342802 : LEES BRONCHUS CLAMP**

Angled, DeBakey Jaw, 1 5/8" (4cm) tip length, Curved shaft, 10" (25.5cm)



**PG342804 : BRONCHUS CLAMP**

Angled, DeBakey Jaw, 3 1/4" (8.3cm) tip length, 9" (23cm)



**PG342822 : SAROT BRONCHUS CLAMP**

8.9cm tip length, 9 1/4" (23.5cm)