

CVT



PG351500 : DEBAKEY SUCTION TUBE

9 mm tip, 11-3/8" (29.0 cm)



PG351501 : SUCTION TUBE TIP FOR 351500

3/8" (1.0 cm)



PG351502 : DEBAKEY SUCTION TUBE TIP

Curved, 3 mm tip, 12-1/4" (31.0 cm)



PG351504 : DEBAKEY-ADSON CORONARY SUCTION TUBE

4.5 mm tip, 6-1/8" (15.4 cm)



PG351505 : MORSE-ANDREWS SUCTION TUBE

6 mm tip, 10-1/4" (26.0 cm)



PG351506 : DEBAKEY-ADSON CORONARY SUCTION TUBE

4 mm tip, 8-1/8" (20.5 cm)

PG351508 : DEBAKEY FIBER OPTIC SUCTION TUBE

Curved, 9 mm diameter tip, 11-3/8" (28.9 cm)



PG351547 : COOLEY SUMP-SUCTION TUBE

Curved, 9 mm tip, 12-7/8" (32.7 cm)



PG351560 : STONEY HEPARIN INJECTOR

1-7/8" (4.5 cm)



PG351573 : ZANGER-DEBAKEY HEPARIN NEEDLE

4 mm tip, 2-1/2" (6.4 cm)



PG351574 : HORSLEY HEPARIN NEEDLE

2-3/8" (5.9 cm)



PG351575 : DEBAKEY HEPARIN NEEDLE SMALL

Malleable, 3 mm tip, 2-1/2" (6.4 cm)



PG351576 : DEBAKEY HEPARIN NEEDLE MEDIUM

Malleable, 3.7 mm tip, 2-7/8" (7.4 cm)



PG351577 : DEBAKEY HEPARIN NEEDLE LARGE

Malleable, 4 mm tip, 6" (15.1 cm)



PG351578 : IMA HEPARIN NEEDLE

23 gauge tip, 2-1/8" (5.3 cm)



PG351586 : GARRETT HEPARIN INJECTOR

Malleable, 5/16" (8 mm) tip, working length 3" (7.5 cm)



PG351587 : GARRETT HEPARIN INJECTOR

Malleable, 1/2" (12 mm) tip, working length 3" (7.5 cm)



PG351593 : SAMSON HEPARIN INJECTOR

Malleable, small 7 mm tip, 3-3/8" (8.5 cm)



PG351594 : SAMSON HEPARIN INJECTOR MEDIUM


Malleable, medium 9 mm tip, 3-3/8" (8.5 cm)



PG351650 : COOLEY VENA CAVA CLAMP


Curved, 28 Fr. ring tip width, 7 1/4" (18.5cm)






PG351651 : COOLEY VENA CAVA CLAMP

Curved, 30 Fr. ring tip width, 7 1/4" (18.5cm)



PG351653 : COOLEY VENA CAVA CLAMP

Curved, 36 Fr. ring tip width, 7 1/4" (18.5cm)



PG351654 : COOLEY VENA CAVA CLAMP

Curved, 34 Fr. ring tip width, 7 1/4" (18.5cm)



PG351664 : MAINI RICHARDSON DOUBLE ENDED CENTER BLADE FOR HENLY RETRACTOR

Blade 31mm & 40mm wide



PG351666 : RETRACTOR SIDE BLADE FOR 351675 PILLING MULTIPURPOSE RETRACTOR

4 Blunt prongs, Side Blade 7/8" (2.2cm) wide x 1 1/8" (2.9cm) deep