



PG351667 : RETRACTOR SIDE BLADE FOR 351675 PILLING MULTIPURPOSE RETRACTOR

4 Blunt prongs, Side Blade 7/8" (2.2cm) wide x 2 1/8" (5.5cm) deep



PG351668 : RETRACTOR SIDE BLADE FOR 351675 PILLING MULTIPURPOSE RETRACTOR

4 Blunt prongs, Side Blade 7/8" (2.2cm) wide x 3 1/8" (8.0cm) deep



PG351669 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 7/8" (2.2cm) deep



PG351670 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 1 1/8" (2.8cm) deep



PG351671 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 1 1/4" (3.2cm) deep



PG351672 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 2" (5.2cm) deep



PG351673 : RETRACTOR CENTER BLADE 351675 PILLING MULTIPURPOSE RETRACTOR

Sharp, Center blade 5/8" (1.6cm) wide x 2 1/2" (6.5cm) deep



PG351675 : PILLING® MULTIPURPOSE RETRACTOR COMPLETE

A versatile, self-retaining Henly-type retractor designed for use in the following surgical specialties: peripheral vascular surgery in the carotid, subclavian, femoral, popliteal and tibial regions. Spinal surgery in the...



PG351676 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 1" (2.5cm) deep, 1 pair



PG351677 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 2" (5cm) deep, 1 pair



PG351678 : IMPROVED HENLY RETRACTOR SIDE BLADE

Side Blade 7/8" (2.2cm) wide x 3" (7.6cm) deep, 1 pair



PG351679 : IMPROVED HENLY RETRACTOR WITH STRAIGHT ARMS

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame



PG351680 : IMPROVED HENLY RETRACTOR WITH ARMS SPREAD OUTWARD

7" (18cm), 3 pair of Side Blades 1" (2.5cm), 2" (5cm) and 3" (7.6cm) deep x 7/8" (2.2cm) wide. 3 Center Blades 7/8" (2.2cm), 1 1/8" (2.8cm), 1 1/4" (3.2cm) deep x 5/8" (1.6cm) wide, frame


PG351715 : CONN AORTIC COMPRESSOR

15 1/2" (39.5cm)




PG351743 : REVUELTA-GARCIA ALLIS CLAMP

Half Curved, DeBakey Jaw, 1" (2.4cm) tip length, 10 1/8" (25.5cm)



PG351744 : REVUELTA-GARCIA ALLIS CLAMP

Half Curved, DeBakey Jaw, 10 7/8" (27.5cm)



PG351745 : REVUELTA-GARCIA VALVE SCISSORS

Curved, 10 7/8" (27.5cm)



PG351747 : REVUELTA GARCIA VALVE LEAFLET RETRACTOR

Curved, 24.0mm width blade



PG351748 : CARDIAC VALVE FORCEPS


Serrated, 3mm tip, 9 1/2" (25 cm)



PG351750 : BLANCO VALVE RETRACTOR

5mm width x 3/8" (1.0cm) depth blade, 11 1/2" (29.0cm)





PG351752 : BLANCO VALVE RETRACTOR

7mm width x 3/8" (1.0cm) depth blade, 11" (28.0cm)




PG351754 : CUSHING-BLANCO VALVE RETRACTOR

12mm width x 1/4" (0.5cm) depth blade, 11" (28.0cm)




PG351755 : BLANCO VALVE RETRACTOR

9mm width x 1/2" (1.0cm) depth blade, 11 1/2" (29.0cm)



PG351763 : ABLAZA-BLANCO RETRACTOR

9/16" (1.4cm) width blade, 13" (33.0cm)



PG351767 : DESMARRE VALVE RETRACTOR

4mm width blade, 7 1/4" (18.5cm)

