

CVT



PG351768 : DESMARRE VALVE RETRACTOR

6mm width blade, 7 1/4" (18.5cm)



PG351769 : DESMARRE VALVE RETRACTOR

8mm width blade, 7 1/4" (18.5cm)



PG351770 : DESMARRE VALVE RETRACTOR

10mm width blade, 7 1/4" (18.5cm)



PG351771 : DESMARRE VALVE RETRACTOR

12mm width blade, 7 1/4" (18.5cm)



PG351772 : DESMARRE VALVE RETRACTOR

15mm width blade, 7 1/4" (18.5cm)



PG351773 : CRAWFORD AORTIC RETRACTOR

15mm width x 3/4" (2.0cm) depth blade, 7 3/8" (18.5cm)




PG351774 : CRAWFORD AORTIC RETRACTOR

22mm width x 3/4" (2.0cm) depth blade, 7 7/8" (20.0cm)



PG351779 : HENLY POPLITEAL RETRACTOR

4 X 3-pronged, Supplied w/3 center blades: 5/8" (1.5cm) width x 3/4" (1.8cm), 5/8" (1.5cm) x 1" (2.5cm) and 5/8" (1.5cm) x 1 1/4" (3.1cm) depth blade, 6 3/8" (16.0cm)




PG351783 : LEMOLE ATRIAL RETRACTOR

1 7/8" (5.0cm) width blade, 10 1/2" (26.5cm)



PG351784 : HENLY DEEP CENTER BLADES RETRACTOR

5/8" (1.6cm) width x 1 1/2" (3.8cm) / 1 3/4" (4.5cm) depth blade, 4 7/8" (12.5cm)



PG351785 : DEBAKEY VALVE HOOKS

Small, Straight, 3/8" (8.6mm) hook, 12 1/2" (31.6cm)



PG351788 : COOLEY MITRAL VALVE RETRACTOR, SMALL

Curved right, 7/8" (2.2cm) width blade, 7 1/4" (18.5cm)



PG351789 : COOLEY MITRAL VALVE RETRACTOR, MEDIUM

Curved right, 1 3/8" (3.4cm) width blade, 7 1/2" (19.0cm)

PG351790 : COOLEY MITRAL VALVE RETRACTOR, LARGE

Curved right, 1 7/8" (4.7cm) width x 1 3/8" (3.5cm) depth blade, 7 3/8" (19.0cm)



PG351791 : COOLEY MITRAL VALVE RETRACTOR

Curved left, 2" (5cm) width blade, 8 7/8" (22.5cm)



PG351800 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 6" (16 cm)



PG351801 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 6" (16 cm)



PG351802 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 7 3/4" (20 cm)



PG351803 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 7 3/4" (20 cm)



PG351803DG : DEBAKEY-DIETHRICH DIAMA-GRIP TISSUE FORCEPS

Straight, 1.3 mm tip, 7-3/4" (20 cm), Diama-Grip jaws





PG351804 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 9 1/2" (24 cm)



PG351806 : DEBAKEY TISSUE FORCEPS

Straight, Delicate, 1mm tip, 9 1/2" (24 cm)



PG351808 : DEBAKEY TISSUE FORCEPS

Straight, 2mm tip, 12" (31 cm)



PG351810 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 6" (16 cm)



PG351812 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 7 3/4" (20 cm)