



PG351814 : DEBAKEY TISSUE FORCEPS

Angled, 2mm tip, 9 1/2" (24 cm)



PG351815 : DEBAKEY TISSUE FORCEPS

2mm tip, 6" (15 cm), TC



PG351817 : DEBAKEY TISSUE FORCEPS

2mm tip, 7 3/4" (20 cm), TC



PG351818 : DEBAKEY TISSUE FORCEPS

Delicate, 7 3/4" (20 cm), TC



PG351819 : DEBAKEY TISSUE FORCEPS

9 1/2" (24 cm), TC



PG351820 : DEBAKEY ULTRALIGHT TISSUE FORCEPS

1.5mm tip, 6" (16 cm)



PG351821 : DEBAKEY ULTRALIGHT TISSUE FORCEPS

1.5mm tip, 7 3/4" (20 cm)



PG351824 : DEBAKEY TISSUE FORCEPS

Delicate, 9 1/2" (24 cm), TC



PG351825 : WHEELER VESSEL FORCEPS

DeBakey type jaw, 1.6mm , 5 1/2" (14 cm)



PG351827 : WHEELER VESSEL FORCEPS

DeBakey type jaw, 7 1/8" (18 cm)




PG351829 : WHEELER VESSEL FORCEPS

DeBakey type jaw, 8 3/4" (23 cm)



PG351830 : GREGORY SUTURE STAY CLAMP

Straight, 3 7/8" (10cm)



PG351831 : GREGORY SUTURE STAY CLAMP

Curved, 3 7/8" (10cm)

PG351835 : ADAMS TISSUE FORCEPS

Delicate, 1mm tip, 8" (21 cm), TC



PG351840 : SELMAN TISSUE FORCEPS

DeBakey type jaw, 6" (15 cm)



PG351841 : SELMAN TISSUE FORCEPS

DeBakey type jaw, 8" (21 cm)



PG351842 : SELMAN TISSUE FORCEPS

DeBakey type jaw, 10" (26 cm)



PG351843 : DEBAKEY BAYONET TISSUE FORCEPS

7 3/4" (20 cm)



PG351845 : DEBAKEY BAYONET TISSUE FORCEPS

Fine, 7 3/4" (20 cm)



PG351850 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 1.3mm tip, 7 1/4" (19 cm)





PG351852 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 8 1/4" (21 cm)



PG351854 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 9 5/8" (25 cm)



PG351856 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 10 1/4" (26 cm)



PG351858 : POTTS SMITH DRESSING FORCEPS

Fine, Serrated, 11 7/8" (30 cm)



PG351860 : CUSHING DRESSING FORCEPS

6 7/8" (18 cm), TC

