

Home // CVT

CVT

NO IMAGE
AVAILABLE**PG353011 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 1 1/2" (4cm)

NO IMAGE
AVAILABLE**PG353012 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 1 7/8" (5cm)

NO IMAGE
AVAILABLE**PG353014 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 2 1/2" (6.5cm)

NO IMAGE
AVAILABLE**PG353016 : JOHNS HOPKINS BULLDOG CLAMP**

Curved, Serrated, 3" (7.5cm)

NO IMAGE
AVAILABLE**PG353025 : HESS BULLDOG CLAMP**

Angled, Serrated, 1 1/8" (2.7cm) tip length, 2 1/8" (5.5cm)

NO IMAGE
AVAILABLE**PG353035 : BULLDOG CLAMP APPLYING FORCEPS**

Angled, Smooth, 9 1/2" (24 cm)

NO IMAGE
AVAILABLE

PG353040 : SCHWARTZ CLAMP HOLDING FORCEPS

Angled, 7 1/2" (19 cm)

NO IMAGE
AVAILABLE

PG353041 : SCHWARTZ MICRO BULLDOG CLAMP

Straight, 1" (2.5cm)

NO IMAGE
AVAILABLE

PG353042 : SCHWARTZ MICRO BULLDOG CLAMP

Slight Curved, 1/4" (6.4mm) tip length, 1" (2.5cm)

NO IMAGE
AVAILABLE

PG353043 : SCHWARTZ MICRO BULLDOG CLAMP

Strongly Curved, 3/8" (10.2mm) tip length, 1" (2.5cm)

NO IMAGE
AVAILABLE

PG353050 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Straight, 3/4" (2cm) tip length, 3 1/8" (8cm)

NO IMAGE
AVAILABLE

PG353052 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Straight, 1 1/8" (3cm) tip length, 3 1/2" (9cm)

NO IMAGE
AVAILABLE

PG353054 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Straight, 2" (5.1cm) tip length, 4 1/4" (11cm)

NO IMAGE
AVAILABLE

PG353056 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Straight, 2 5/8" (6.8cm) tip length, 4 7/8" (12.5cm)

NO IMAGE
AVAILABLE

PG353060 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 5/8" (1.6cm) tip length, 2 7/8" (7.5cm)

NO IMAGE
AVAILABLE

PG353062 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 1" (2.5cm) tip length, 3 1/4" (8cm)

NO IMAGE
AVAILABLE

PG353064 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 1 5/8" (4.1cm) tip length, 3 7/8" (10cm)

NO IMAGE
AVAILABLE

PG353066 : DEBAKEY CROSS ACTION BULLDOG CLAMP

Curved, 2 1/8" (5.3cm) tip length, 4 1/4" (11cm)

NO IMAGE
AVAILABLE

PG353070 : SANTULLI BULLDOG CLAMP

Angled, DeBaKey Jaw, 1 5/8" (4.1cm) tip length, 2 3/4" (7cm)

NO IMAGE
AVAILABLE

PG353075 : WYLIE LUMBAR BULLDOG CLAMP

Slight Curved, DeBaKey Jaw, 7/8" (2.1cm) tip length, 3" (7.5cm)

AVAILABLE

NO IMAGE
AVAILABLE

PG353077 : GREGORY CAROTID BULLDOG CLAMP

Angled, DeBakey Jaw, 5/8" tip length, 4 1/2" (11.5cm)

NO IMAGE
AVAILABLE

PG353078 : GREGORY CAROTID BULLDOG CLAMP

Angled, DeBakey Jaw, 1/2" tip length, 4 1/2" (11.5cm)

NO IMAGE
AVAILABLE

PG353080 : CALMAN CAROTID BULLDOG CLAMP

Straight, DeBakey Jaw, 1" (2.4cm) tip length

NO IMAGE
AVAILABLE

PG353085 : CALMAN CAROTID BULLDOG CLAMP

Angled, DeBakey Jaw, 1 5/8" (4cm) tip length, 3 3/4" (9.5cm)

NO IMAGE
AVAILABLE

PG353086 : JAVID CAROTID SHUNT CLAMP

Curved, Longitudinal Serrations, 1 3/4" (4.45cm) tip length, 7 3/8" (18.5cm)