




PG353155 : DIETHRICH BULLDOG CLAMP

Angled, Serrated, 3/4" (20mm) tip length, 2 1/8" (5.5cm)



PG353182 : DEBAKEY RING HANDLE BULLDOG CLAMP

Straight, 2 1/4" (5.7cm) tip length, 2.4mm tip width, Straight shaft, 5 3/8" (13.5cm)



PG353183 : GREGORY RING HANDLE BULLDOG CLAMP

Straight, DeBakey Jaw, 1.6mm tip width, Straight shaft, 4 3/4" (12cm)



PG353186 : DEBAKEY RING HANDLE BULLDOG CLAMP

Half Curved, 2 3/8" (5.9cm) tip length, 2.4mm tip width, Curved shaft, 5 3/8" (13.5cm)



PG353187 : GREGORY RING HANDLE BULLDOG CLAMP

Half Curved, DeBakey Jaw, 1 1/8" (2.7cm) tip length, 1.5mm tip width, Curved shaft, 4 1/8" (10.5cm)



PG353189 : DEBAKEY RING HANDLE BULLDOG CLAMP

Angled 45, 2 3/8" (6cm) tip length, 2.4mm tip width, Curved shaft, 5 1/4" (13.5cm)




PG353190 : GREGORY RING HANDLE BULLDOG CLAMP

Angled 45, DeBakey Jaw, 1 5/8" (4.2cm) tip length, 1.5mm tip width, Curved shaft, 4 3/8" (11cm)




PG353192 : DEBAKEY RING HANDLE BULLDOG CLAMP

Angled 90, 1 3/4" (4.6cm) tip length, 2.4mm tip width, Straight shaft, 4 1/8" (10.5cm)




PG353195 : DEBAKEY RING HANDLE BULLDOG CLAMP

Slight Curved, 2 1/2" (6.3cm) tip length, 2.4mm tip width, Curved shaft, 5 1/8" (13cm)




PG353203 : STYLET FOR RUMEL TOURNIQUET

10 1/4" (26cm)




PG353206 : STYLET FOR WECK STYLE RUMEL TOURNIQUET

9 3/8" (24cm)



PG353207 : STYLET CHILD SIZE RUMEL TOURNIQUET

7" (18cm)



PG353490 : LELAND JONES CLAMP

Straight, DeBakey Jaw, 3 3/8" (8.7cm) tip length, 2.4mm tip width, 8 1/8" (20.5cm)

PG353491 : LELAND JONES CLAMP

Angled 15, DeBakey Jaw, 3 1/2" (8.9cm) tip length, 2.4mm tip width, 8 1/8" (20.5cm)

PG353492 : LELAND JONES CLAMP

Angled 30, DeBakey Jaw, 3 5/8" (9.05cm) tip length, 2.4mm tip width, 8" (20.5cm)

PG353494 : GREGORY PROFUNDA CLAMP

DeBakey Jaw, 2 1/4" (5.8cm) tip length, 1.7mm tip width, 4 3/4" (12cm)

PG353495 : GREGORY PROFUNDA CLAMP

DeBakey Jaw, 3 1/4" (8.3cm) tip length, 2.8mm tip width, 6 1/8" (15.5cm)

PG353496 : GREGORY PROFUNDA CLAMP



DeBakey Jaw, 4" (10.2cm) tip length, 2.8mm tip width, 7" (18cm)

PG353500 : DEBAKEY COARCTATION CLAMP

Straight, 3 1/8" (7.9cm) tip length, 2.8mm tip width, 9" (23cm)


PG353503 : DEBAKEY COARCTATION CLAMP

Straight, 3" (7.6cm) tip length, 2.8mm tip width, Angled shaft, 8 3/4" (22cm)




PG353510 : DEBAKEY PATENT DUCTUS CLAMP

Straight, 2 1/4" (5.7cm) tip length, 2.8mm tip width, Straight shaft, 8 1/4" (21cm)




PG353513 : DEBAKEY PATENT DUCTUS CLAMP

Straight, 2 3/8" (5.9cm) tip length, 2.8mm tip width, Angled shaft, 8 1/8" (20.5cm)




PG353515 : DEBAKEY RENAL ARTERY CLAMP

3 1/2" 9cm tip length, Curved shaft, 10 1/8" (25.5cm)



PG353522 : WIRTHLIN SPLENO RENAL SHUNT CLAMP

3 1/8" 7.9cm tip length, 1.7mm tip width, Curved shaft, 6 1/2" (16.5cm)



PG353524 : STONEY SPLENO RENAL SHUNT CLAMP

DeBakey Jaw, 3 1/4" (8.3cm) tip length, 2mm tip width, Curved shaft, 9 1/4" (23.5cm)