




PG353528 : DALE FEM-POP CLAMP

DeBakey Jaw, 4 1/8" (10.4cm) tip length, 2.8mm tip width, Curved shaft, 7 3/8" (19cm)



PG353531 : CLEVELAND CLINIC RENAL CLAMP

Left Curved, DeBakey Jaw, 3 3/8" (8.6cm) tip length, 9 5/8" (24.5cm)



PG353532 : CLEVELAND CLINIC RENAL CLAMP

Right Curved, DeBakey Jaw, 3 3/8" (8.6cm) tip length, 9 5/8" (24.5cm)



PG353533 : CLEVELAND CLINIC BABY RENAL CLAMP

Left Curved, DeBakey Jaw, 2mm tip width, 5 7/8" (15cm)



PG353534 : CLEVELAND CLINIC BABY RENAL CLAMP

Right Curved, DeBakey Jaw, 2mm tip width, 5 7/8" (15cm)




PG353535 : DEBAKEY PERIPHERAL VASCULAR CLAMP

ANGLED, 3 3/4" (9.4cm) tip length, 2.8mm tip width, Curved shaft, 7 1/8" (18cm)




PG353536 : DIETHRICH CAVAL OCCLUDING CLAMP

Curved, DeBakey Jaw, 4 1/4" (10.7cm) tip length, Curved shaft, 9 3/4" (25cm)




PG353538 : HENLY SUBCLAVIAN ARTERY CLAMP

4 1/8" 10.4cm tip length, 2.8mm tip width, 8 3/8" (21cm)




PG353539 : CHERRY SUPRACELIAC CLAMP

Curved, 2 x 3 AtrauGrip Jaw, Curved shaft, 12 1/4" (31cm)



PG353540 : DEBAKEY "S" PV CLAMP

Curved, 4 1/8" (10.47cm) tip length, 2.8mm tip width, Curved shaft, 7 7/8" (20cm)




PG353541 : DALE FEM-POP CLAMP

Curved, DeBakey Jaw, 3" (7.5cm) tip length, 2.4mm tip width, Curved shaft, 7" (17.5cm)



PG353542 : DALE FEM-POP CLAMP

Curved, DeBakey Jaw, 3 3/8" (8.5cm) tip length, 2.4mm tip width, Curved shaft, 7 1/8" (18cm)




PG353543 : DALE FEM-POP CLAMP

Curved, DeBakey Jaw, 3 1/2" (8.9cm) tip length, 2.4mm tip width, Curved shaft, 7 1/4" (18.5cm)

PG353545 : DEBAKEY AORTIC CLAMP

Curved, 3 7/8" (9.8cm) tip length, 2.4mm tip width, Curved shaft, 7 1/8" (18cm)




PG353546 : FAVALORO PROXIMAL ANASTOMOSIS CLAMP

Curved, DeBakey Jaw, 4 1/8" (10.4cm) tip length, 2.8mm tip width, 8" (20cm)



PG353550 : DEBAKEY PV CLAMP

Curved, 2.3mm tip width, 6 3/4" (17cm)



PG353560 : DEBAKEY BAINBRIDGE CLAMP

Straight, 2.4mm tip width, 7" (17.5cm)



PG353562 : DEBAKEY BAINBRIDGE CLAMP

Straight, 3 1/8" (8cm) tip length, 2.8mm tip width, 7 7/8" (20cm)



PG353565 : DEBAKEY BAINBRIDGE CLAMP


Curved, 2 3/8" (5.9cm) tip length, 2.4mm tip width, 7" (17.5cm)



PG353567 : DEBAKEY BAINBRIDGE CLAMP


Curved, 3 1/4" (8.1cm) tip length, 2.4mm tip width, 7 5/8" (19.5cm)





PG353568 : DEMOS TIBIAL ARTERY CLAMP

Straight, DeBakey Jaw, 2 3/8" (6cm) tip length, 2.1mm tip width, Angled shaft, 6 3/4" (17cm)




PG353569 : WYLIE CLOSER CAROTID CLAMP

DeBakey Jaw, 2 1/2" (6.3cm) tip length, 1.7mm tip width, Curved shaft, 6 1/8" (15.5cm)



PG353570 : WYLIE EXTERNAL CAROTID CLAMP

DeBakey Jaw, 3 3/4" (9.5cm) tip length, 2.8mm tip width, Curved shaft, 8" (20.5cm)



PG353571 : WYLIE INTERNAL CAROTID CLAMP

DeBakey Jaw, 2.7mm tip width, Curved shaft, 8 1/4" (21cm)



PG353580 : KAY ANASTOMOSIS AORTIC CLAMP

Curved, DeBakey Jaw, 1 1/4" (3.2cm) tip length, Angled shaft, 7 5/8" (19.5cm)