



PG354730 : DEBAKEY PEDIATRIC CLAMP

Curved, 1 3/4" (4.4cm) tip length, 2.8mm tip width, Curved shaft, 6 7/8" (17.5cm)



PG354735 : DEBAKEY MULTIPURPOSE PEDIATRIC CLAMP

Angled 30, 3/4" (1.9cm) tip length, 2.8mm tip width, Curved shaft, 6 3/4" (17cm)




PG354737 : DEBAKEY MULTIPURPOSE PEDIATRIC CLAMP

Angled 60, 1 1/4" (3.2cm) tip length, 2.7mm tip width, Curved shaft, 6 1/4" (16cm)




PG354739 : DEBAKEY MULTIPURPOSE PEDIATRIC CLAMP

Angled 90, 1 1/2" (3.7cm) tip length, 2.8mm tip width, Curved shaft, 6" (15cm)



PG354750 : BECK MINIATURE PEDIATRIC CLAMP

Angled, DeBakey Jaw, 1 1/8" (2.86cm) tip length, 2.8mm tip width, Straight shaft, 7" (18cm)




PG354755 : GLOVER MINIATURE PEDIATRIC CLAMP

Angled, DeBakey Jaw, 1/2" (1.3cm) tip length, 2.8mm tip width, Straight shaft, 7 1/8" (18cm)



PG354800 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Straight, 1 7/8" (4.9cm) tip length, 2mm tip width, Straight shaft, 6 3/4" (17cm)



PG354801 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Straight, 1 7/8" (4.9cm) tip length, 2mm tip width, Angled shaft, 6 1/2" (16.5cm)



PG354802 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Straight, 2 1/4" (5.8cm) tip length, 2mm tip width, Straight shaft, 7" (18cm)



PG354803 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Straight, 2 3/8" (6cm) tip length, 2mm tip width, Angled shaft, 7" (17.5cm)




PG354804 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Slight Curved, 2 5/8" (6.8cm) tip length, 2mm tip width, Straight shaft, 7 1/4" (18.5cm)



PG354805 : COOLEY PEDIATRIC MULTIPURPOSE CLAMP

Angled 30, 3/4" (1.9cm) tip length, 2mm tip width, Angled shaft, 6 3/4" (17cm)




PG354806 : COOLEY PEDIATRIC MULTIPURPOSE CLAMP

Angled 60, 1 1/4" (3.3cm) tip length, 2mm tip width, Curved shaft, 6 1/2" (16.5cm)

PG354807 : COOLEY PEDIATRIC MULTIPURPOSE CLAMP

Angled 90, 1 1/2" (3.8cm) tip length, 2mm tip width, Curved shaft, 5 3/4" (14.5cm)




PG354808 : COOLEY PEDIATRIC PATENT DUCTUS CLAMP

Spoon, 2 1/2" (6.2cm) tip length, 2mm tip width, Straight shaft, 7 1/8" (18cm)



PG354810 : COOLEY PEDIATRIC CLAMP

Curved, 1" (2.5cm) tip length, 2mm tip width, Straight shaft, 6 1/8" (15.5cm)



PG354815 : COOLEY PEDIATRIC CLAMP

Curved, 3/4" (1.9cm) tip length, 2mm tip width, Curved shaft, 7" (18cm)



PG354816 : COOLEY PEDIATRIC CLAMP

Curved, 1 1/8" (2.7cm) tip length, 2mm tip width, Curved shaft, 7" (18cm)



PG354817 : COOLEY PEDIATRIC CLAMP

Curved, 1 1/8" (3cm) tip length, 2mm tip width, Curved shaft, 7 1/8" (18cm)



PG354826 : COOLEY PEDIATRIC CLAMP

Full Curved, 1 1/2" (3.8cm) tip length, 2mm tip width, Angled shaft, 6 5/8" (17cm)





PG354827 : COOLEY PEDIATRIC CLAMP

Strongly Curved, 1 3/4" (4.6cm) tip length, 2mm tip width, Curved shaft, 6 3/4" (17cm)




PG354830 : COOLEY PEDIATRIC CLAMP

Angled 90, small jaws, 6 3/4" (17cm)




PG354835 : COOLEY PEDIATRIC CLAMP

Angular shafts, small double angled 60 calibrated jaw 1 1/8" (2.8cm), spring handle, 1 x 7 ratchets, 7 1/8" (18cm)



PG354836 : COOLEY PEDIATRIC CLAMP

Angular shafts, large double angled 60 calibrated jaw 1 1/4" (3.2cm), spring handle, 1 x 7 ratchets, 7 1/8" (18cm)



PG354839 : COOLEY RENAL CLAMP

Angled, 5/8" (1.7cm) tip length, 2mm tip width, Curved shaft, 7 3/4" (19.5cm)

