

**PG354845 : COOLEY ANASTOMOSIS CLAMP**

Angled, 3/4" (2cm) tip length, 2mm tip width, 7 1/8" (18cm)



**PG354847 : COOLEY ANASTOMOSIS CLAMP**

Angled 30, 5/8" (1.6cm) tip length, 2mm tip width, 7 1/8" (18cm)



**PG354900 : DEBAKEY NEEDLE HOLDER**

13.7 x 2 mm jaw, 7" (18 cm), TC




**PG354902 : DEBAKEY NEEDLE HOLDER**

17 x 1.8 mm jaw, 9-1/4" (23.5 cm), TC




**PG354904 : DEBAKEY NEEDLE HOLDER**

21.6 x 2.3mm jaw, 10 1/2" (26.5cm), TC




**PG354906 : DEBAKEY NEEDLE HOLDER**

19 x 2.6mm jaw, 11 3/4" (30cm), TC




**PG354915 : W.M.LEMMON NEEDLE HOLDER**

1.5mm jaw, 5 3/4" (14.5cm), TC




**PG354916 : W.M.LEMMON NEEDLE HOLDER**

1.5mm jaw, 6 3/4" (17cm), TC




**PG354917 : W.M. LEMMON NEEDLE HOLDER**

7.2 x 1.6mm jaw, 7 1/2" (19cm), TC




**PG354920 : HOWELL CORONARY NEEDLE HOLDER**

14 x 1.6mm jaw, 7 1/4" (18.5cm), TC



**PG354928 : HOWELL VALVE NEEDLE HOLDER**

15.3 x 2mm jaw, 8" (20.5cm), TC



**PG354931 : BAILEY NEEDLE HOLDER**

17 x 2mm jaw, 7" (18.5cm), TC



**PG354941 : JAMISON NEEDLE HOLDER**

1.8mm jaw, 9" (23cm), TC



**PG354958 : PILLING THIN VASCULAR NEEDLE HOLDER**

16.3 x 1.8mm jaw, 9" (23cm), TC



**PG354959 : PILLING THIN VASCULAR NEEDLE HOLDER**

16 x 2.2mm jaw, 10 1/4" (26cm), TC



**PG354960 : PILLING THIN VASCULAR NEEDLE HOLDER**

11 x 1.8mm jaw, 6" (15cm), TC



**PG354961 : PILLING THIN VASCULAR NEEDLE HOLDER**

11.8 x 1.7mm jaw, 7" (17.5cm), TC



**PG354962 : PILLING THIN VASCULAR NEEDLE HOLDER**

15.8 x 1.8 mm jaw, 8-1/4" (21 cm), TC



**PG354964 : CASTROVIEJO TECHNIQUE® NEEDLE HOLDERS**

Straight, with lock, 13 x 1.4 mm serrated jaw, 7-1/4" (18.5 cm), TC



**PG354964DG : CASTROVIEJO NEEDLE HOLDER**

Straight, Flat handle, With lock, 1.8mm Diamo-Grip jaw, 7 1/4" (18.5cm)



---

### **PG354965 : CASTROVIEJO TECHNIQUE® NEEDLE HOLDERS**

Straight, with lock, 13 x 1.4 mm serrated jaw, 8" (20.5 cm), TC



---

### **PG354965DG : CASTROVIEJO NEEDLE HOLDER**

Flat handle, w/Lock, Diamo-Grip, 8" (20.5cm)



---

### **PG354966 : MICRO SURGERY TECHNIQUE® NEEDLE HOLDERS**

Straight, Flat handle, With lock, 13 x 1.4mm smooth jaw, 7" (18cm), TC



---

### **PG354967 : MICRO SURGERY TECHNIQUE® NEEDLE HOLDERS**

Straight, flat handle, with lock, 13 x 1.4 mm smooth jaw, 8" (20.5 cm), TC



---

### **PG354969 : SPENCER CORONARY TECHNIQUE® NEEDLE HOLDER**

Straight, smooth, 12 x 1.2 mm jaw, 0.9 mm tip, 7" (18 cm), TC

