



PG354971 : GIANNINI NEEDLE HOLDER

Straight, with lock, 1.8 mm jaw, 10-1/4" (26 cm), TC



PG354972 : CASTROVIEJO NEEDLE HOLDER

Straight, Round handle, With lock, 11.5 x 0.8mm Diamo-Grip jaw, 6" (15cm)



PG354973 : CASTROVIEJO NEEDLE HOLDER

Curved, Round handle, With lock, 0.9mm Diamo-Grip jaw, 6" (15cm)




PG354974 : CASTROVIEJO NEEDLE HOLDER

Straight, Round handle, With lock, 0.8mm Diamo-Grip jaw, 7 1/4" (18.5cm)



PG354975 : CASTROVIEJO NEEDLE HOLDER

Straight, round handle, with lock, 0.7 mm Diamo-Grip jaw, 8-1/4" (21 cm)




PG354981 : PILLING EXTRA LIGHT NEEDLE HOLDER

12.3 x 1.8mm jaw, 6" (15cm), TC




PG354982 : PILLING EXTRA LIGHT NEEDLE HOLDER

14.5 x 1.6 mm jaw, 7" (18 cm), TC




PG354985 : PILLING INTRACARDIAC NEEDLE HOLDER

17 x 1.8mm jaw, 7" (18cm), TC




PG354986 : PILLING INTRACARDIAC NEEDLE HOLDER

14.5 x 1.6mm jaw, 8" (20.5cm), TC



PG354987 : PILLING INTRACARDIAC NEEDLE HOLDER

16 x 1.9mm jaw, 9 1/4" (23.5cm), TC



PG354988 : PILLING INTRACARDIAC NEEDLE HOLDER

16.5 x 2mm jaw, 10 1/4" (26cm), TC



PG354990 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)

11.5 x 1.3mm jaw, 5 1/4" (13.5cm), TC



PG354992 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)

12.7 x 1.2mm jaw, 6" (15.5 cm), TC


PG354994 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)

14.3 x 1.2 mm jaw, 7" (18 cm), TC



PG354996 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)

14.5 x 1.2mm jaw, 8" (20.5cm), TC



PG354998 : RYDER EXTRA FINE NEEDLE HOLDER (FRENCH EYE)

15 x 1.3mm jaw, 9" (23cm), TC



PG355080 : GARRETT VASCULAR DILATOR

Short 1.0 mm, malleable shaft, 5 1/2" (14cm)



PG355081 : GARRETT VASCULAR DILATOR

Short 1.5 mm, malleable shaft, 5 1/2" (14cm)



PG355082 : GARRETT VASCULAR DILATOR

Short 2.0 mm, malleable shaft, 5 1/2" (14cm)



PG355083 : GARRETT VASCULAR DILATOR

Short 2.5 mm, malleable shaft, 5 1/2" (14cm)



PG355084 : GARRETT VASCULAR DILATOR

Short 3.0 mm, malleable shaft, 5 1/2" (14cm)

PG355085 : GARRETT VASCULAR DILATOR

Short 3.5 mm, malleable shaft, 5 1/2" (14cm)

PG355086 : GARRETT VASCULAR DILATOR

Short 4.0 mm, malleable shaft, 5 1/2" (14cm)

PG355087 : GARRETT VASCULAR DILATOR

Short 4.5 mm, malleable shaft, 5 1/2" (14cm)

PG355088 : GARRETT VASCULAR DILATOR

Short 5.0 mm, malleable shaft, 5 1/2" (14cm)