



PG355089 : GARRETT VASCULAR DILATORS 9/SET

1, 1.5, 2, 2.5, 3, 3.5, 4, 4.5 and 5mm



PG355090 : GARRETT VASCULAR DILATORS 9/SET

1, 1.5, 2, 2.5, 3, 3.5, 4, 4.5 and 5mm



PG355091 : GARRETT VASCULAR DILATOR

Long 1.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355092 : GARRETT VASCULAR DILATOR

Long 1.5 mm, malleable shaft, 8-1/4" (21 cm)



PG355093 : GARRETT VASCULAR DILATOR

Long 2.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355094 : GARRETT VASCULAR DILATOR

Long 2.5 mm, malleable shaft, 8-1/4" (21 cm)




PG355095 : GARRETT VASCULAR DILATOR

Long 3.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355096 : GARRETT VASCULAR DILATOR

Long 3.5 mm, malleable shaft, 8-1/4" (21 cm)



PG355097 : GARRETT VASCULAR DILATOR

Long 4.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355098 : GARRETT VASCULAR DILATOR

Long 4.5 mm, malleable shaft, 8-1/4" (21 cm)



PG355099 : GARRETT VASCULAR DILATOR

Long 5.0 mm, malleable shaft, 8-1/4" (21 cm)



PG355100 : GARRETT DILATOR HOLDER

8-1/8" (20.5 cm)



PG355101 : GARRETT DILATOR HOLDER

5-1/8" (13 cm)

PG355175 : COARCTATION HOOK

Blunt, Right-angled, 7mm hook, 10 1/8" (25.7cm)

PG355180 : OLIVEKRONA DOUBLE END DISSECTOR

Light weight, 7 1/4" (18.5cm)

PG355182 : OLIVEKRONA DOUBLE END DISSECTOR

Light weight, 9 1/2" (24cm)

PG355184 : OLIVEKRONA DOUBLE END DISSECTOR

Heavy weight, 7 3/4" (19.5cm)

PG355186 : OLIVEKRONA DOUBLE END DISSECTOR

Heavy weight, 9 1/2" (24cm)

PG355191 : LEMMON INTIMA DISSECTOR

Slight Curve, 6 3/8" (16.5cm)

PG355192 : LEMMON INTIMA DISSECTOR

Acute Curve, 6 3/8" (16.5cm)





PG355196 : WYLIE CAROTID CLAMP

Curved, Cross-serrated, 5/8" (1.55cm) tip length, 2mm tip width, Curved shaft, 5 7/8" (15cm)



PG355275 : HIEBERT DILATOR

0.5mm x 1mm, Malleable shafts, 6" (15.5cm)



PG355276 : HIEBERT DILATOR

1.5mm x 2mm, Malleable shafts, 6" (15.5cm)



PG355277 : HIEBERT DILATOR

2.5mm x 3mm, Malleable shafts, 6" (15.5cm)



PG355278 : HIEBERT DILATOR

3.5mm x 4mm, Malleable shafts, 6" (15.5cm)

