



CVT




PG355490 : BULLET TIP FOR 355492 GARRETT VIEN PASSER

1 1/4" (3cm)



PG355492 : GARRETT VEIN PASSER

22 1/8" (56.5cm)



PG355493 : HALLMAN TUNNELING FORCEPS

Curved, Serrated, 14 5/8" (37 cm)



PG355495 : DEBAKEY TUNNELER

20 3/4" (52.5cm)



PG355497 : CRAWFORD-COOLEY TUNNELER

17 1/2" (44.5cm)



PG355500 : BERRY ANGIO ACCESS TUNNELER

10 3/4" (27.5cm)

NO IMAGE
AVAILABLE

PG355558DG : CASTROVIEJO MICRO NEEDLE HOLDER

Round handle, w/Lock, 1mm Diamo-Grip jaw, 9" (23cm)



PG355559DG : CASTROVIEJO MICRO NEEDLE HOLDER

Round handle, w/Lock, 1mm Diamo-Grip jaw, 10" (25.5cm)

NO IMAGE
AVAILABLE

PG355560DG : CASTROVIEJO MICRO NEEDLE HOLDER

Round handle, w/Lock, Diamo-Grip, 12" (30.5cm)



PG355798 : DIETRICH CLAMP

Right Angle, Serrated, 3/8" (0.9cm) tip length, 1.2mm tip width, 5 3/8" (13.5cm)



PG355799 : DIETRICH CLAMP

Right Angle, Serrated, 1/4" (0.77cm) tip length, 1.3mm tip width, 4 7/8" (12.5cm)



PG355800 : DIETRICH CLAMP

Right Angle, Serrated, 3/8" (1.1cm) tip length, 1.1mm tip width, 7 1/4" (18.5cm)



PG355801 : DIETRICH CLAMP

Right Angle, Longitudinal Serrations, 3/8" (1.1cm) tip length, 1.8mm tip width, 7 1/4" (18.5cm)

PG355810 : EFFLER GROVES DOUBLE END DISSECTOR AND HOOK

Light weight, 12 1/4" (31cm)

PG355819 : EFFLER-GROVES BULLDOG CLAMP

Straight, DeBakey Jaw, 1 1/4" (3.2cm) tip length, 1.8mm tip width, Straight shaft, 2 3/4" (7cm)

PG356000 : COOLEY TISSUE FORCEPS

6 1/4" (16 cm)

PG356002 : COOLEY TISSUE FORCEPS

8" (21 cm)

PG356003 : COOLEY TISSUE FORCEPS

9 1/2" (24 cm)

PG356025 : COOLEY PERIPHERAL VAS CLAMP

Angled, 1 7/8" (4.9cm) tip length, 2.7mm tip width, Curved shaft, 7 1/2" (19cm)

PG356027 : COOLEY ILIAC CLAMP


Angled, 2 1/4" (5.7cm) tip length, Curved shaft, 8 1/4" (21cm)





PG356028 : COOLEY ILIAC CLAMP

Straight, 2 3/4" (7cm) tip length, Angled shaft, 9 3/4" (25cm)



PG356031 : COOLEY PERIPHERAL VAS CLAMP

Curved, 1 1/8" (2.9cm) tip length, 2mm tip width, Straight shaft, 6 3/4" (17cm)



PG356040 : COOLEY PAT-DUCTUS CLAMP

Straight, 2 1/8" (5.3cm) tip length, 2.3mm tip width, Straight shaft, 8" (20.5cm)



PG356041 : COOLEY PAT-DUCTUS CLAMP

Straight, 1 7/8" (4.9cm) tip length, 2.3mm tip width, Angled shaft, 7 7/8" (20cm)



PG356045 : COOLEY COARCTATION CLAMP

Straight, 2 7/8" (7.4cm) tip length, 2.4mm tip width, Straight shaft, 8 3/4" (22.5cm)

