




CVT

**PG356046 : COOLEY COARCTATION CLAMP**


Straight, 2 3/8" (6cm) tip length, 2.4mm tip width, Angled shaft, 8 7/8" (22.5cm)

**PG356051 : COOLEY PARTIAL-OCCL CLAMP**

Angled, 1" (2.6cm) tip length, 2.3mm tip width, Curved shaft, 8 1/8" (20.5cm)

**PG356052 : COOLEY PARTIAL-OCCL CLAMP**


Angled, 1 1/2" (3.9cm) tip length, 2.4mm tip width, Curved shaft, 8 7/8" (22.5cm)

**PG356053 : COOLEY PARTIAL-OCCL CLAMP**

Angled, 1 7/8" (4.8cm) tip length, 2.4mm tip width, Curved shaft, 8 7/8" (22.5cm)

**PG356060 : COOLEY AORTA CLAMP**

Strongly Curved, 2" (5.1cm) tip length, Straight shaft, 11 1/4" (28.5cm)

**PG356070 : COOLEY MULTIPURPOSE CLAMP**

Angled, 1 1/4" (3.3cm) tip length, 2.3mm tip width, Curved shaft, 9 1/8" (23cm)



PG356071 : COOLEY MULTIPURPOSE CLAMP

Angled, 1 5/8" (4.1cm) tip length, 2.3mm tip width, Curved shaft, 8 1/2" (21.5cm)



PG356072 : COOLEY MULTIPURPOSE CLAMP

Angled, 2 1/8" (5.4cm) tip length, 2.3mm tip width, Curved shaft, 8 1/2" (21.5cm)



PG357100 : PEDIATRIC CLAMP

Straight, Serrated, 1 1/2" (3.8cm) tip length, 1.1mm tip width, Straight shaft, 4 5/8" (12cm)



PG357102 : PEDIATRIC CLAMP

Angled, Serrated, 1 1/2" (3.8cm) tip length, 1.1mm tip width, Straight shaft, 4 5/8" (11.5cm)



PG357106 : PEDIATRIC CLAMP

Angled, Serrated, 1 1/8" (2.8cm) tip length, 1.1mm tip width, Straight shaft, 4" (10.5cm)



PG357110 : PEDIATRIC CLAMP

Slight Curved, 1/2" (1.18cm) tip length, 1.1mm tip width, Straight shaft, 4 1/2" (11.5cm)




PG357180 : CALMAN PEDIATRIC CLAMP

Angled, Serrated, 1" (2.6cm) tip length, 1.1mm tip width, Curved shaft, 4 1/2" (11.5cm)



PG357185 : EXAMPLE PEDIATRIC CLAMP

Angled, Serrated, 1" (2.5cm) tip length, 1.1mm tip width, Curved shaft, 5 5/8" (14cm)



PG357191 : CASTANEDA PEDIATRIC CLAMP

Curved, Serrated, 1/2" (1.4cm) tip length, 1.1mm tip width, Curved shaft, 6" (15cm)




PG357192 : CASTANEDA PEDIATRIC CLAMP

Curved, Serrated, 3/4" (2.03cm) tip length, 1.1mm tip width, Curved shaft, 6" (15.5cm)



PG357194 : THOMPSON PEDIATRIC CLAMP

Angled, Serrated, 1/2" (1.15cm) tip length, 1.1mm tip width, Curved shaft, 6 1/8" (15.5cm)



PG357195 : THOMPSON CAROTID CLAMP

Angled, DeBakey Jaw, 1.6mm tip width, Curved shaft, 5 3/4" (14.5cm)



PG357197 : KARTCHNER CAROTID CLAMP

Angled, Serrated, 3/8" (1.1cm) tip length, 1.1mm tip width, Curved shaft, 6 3/8" (16.5cm)




PG357198 : KARTCHNER CAROTID CLAMP

Angled, DeBakey Jaw, 1.7mm tip width, Curved shaft, 6 1/4" (16cm)



PG357200 : PEDIATRIC CLAMP

Straight, Serrated, 1 7/8" (4.8cm) tip length, 1.1mm tip width, Straight shaft, 6 1/8" (15.5cm)



PG357202 : PEDIATRIC CLAMP

Angled, Serrated, 1/2" (1.27cm) tip length, 1.1mm tip width, Straight shaft, 6" (15.5cm)



PG357262 : PEDIATRIC CLAMP

Angled, Serrated, 3/4" (1.9cm) tip length, 1.1mm tip width, Curved shaft, 5 5/8" (14.5cm)



PG357270 : PEDIATRIC CLAMP

Angled, Serrated, 3/4" (1.8cm) tip length, 1.1mm tip width, Curved shaft, 5 7/8" (15cm)



PG357302 : COOLEY NEONATAL CLAMP

Straight, 2mm tip width, Angled shaft, 5 1/8" (13cm)