



PG357304 : COOLEY NEONATAL CLAMP

Angled, 2mm tip width, Curved shaft, 5" (12.5cm)



PG357306 : COOLEY NEONATAL CLAMP

Curved, 3/4" (19mm) tip length, 2mm tip width, Straight shaft, 4 5/8" (11.5cm)



PG357308 : COOLEY NEONATAL CLAMP

Angled, 2mm tip width, Curved shaft, 4 5/8" (12cm)



PG357310 : COOLEY NEONATAL CLAMP

Angled, 2mm tip width, Curved shaft, 4 3/4" (12cm)



PG357532 : MICRO VASCULAR NERVE HOOK

Small, 7 1/4" (18.4cm)



PG357580 : MICRO BULLDOG CLAMP

Straight, Cross-serrated, 1/4" (7mm) tip length, 1mm tip width, 3/4" (2cm)



PG357581 : MICRO BULLDOG CLAMP

Straight, Cross-serrated, 3/8" (8mm) tip length, 1.3mm tip width, 7/8" (2cm)



PG357582 : MICRO BULLDOG CLAMP

Angled, Cross-serrated, 3/8" (9mm) tip length, 1.4mm tip width, 7/8" (2cm)



PG357583 : MICRO BULLDOG CLAMP

Curved, Cross-serrated, 3/8" (9.5mm) tip length, 1.2mm tip width, 3/4" (2cm)



PG357586 : MICRO BULLDOG BAYONET CLIP APPLYING FORCEPS

7 1/2" (19 cm)



PG357590 : MINI BULLDOG CLAMP

Straight, Cross-serrated, 1/8" (4mm) tip length, 1mm tip width, 1/2" (1cm)



PG357591 : MINI BULLDOG CLAMP

Straight, Cross-serrated, 1/4" (6.3mm) tip length, 1mm tip width, 5/8" (1.5cm)



PG357592 : MINI BULLDOG CLAMP

Curved, Cross-serrated, 1/4" (7mm) tip length, 5/8" (1.5cm)


PG357593 : MINI-ANEURYSM BAYONET CLIP APPLYING FORCEPS

8 1/4" (21 cm)



PG357594 : YASARGIL VASCULAR SCISSORS

Angled 45, bayonet, 7-1/2" (19 cm), 15 mm blades



PG357595 : YASARGIL VASCULAR SCISSORS

Reverse angled 125, Bayonet, 7 1/4" (18.5cm), 12.0mm blades



PG357601 : JACOBSON NEEDLE HOLDER

Straight, Flat handle, With lock, 13 x 1.3mm smooth jaw, 7" (18cm)



PG357603 : JACOBSON NEEDLE HOLDER

Straight, Flat handle, With lock, 13.5 x 1.2mm smooth jaw, 8 1/2" (21.5cm)



PG357603DG : JACOBSON MICRO NEEDLE HOLDER

Flat handle, w/Lock, 18 x 1.3mm Diama-Grip jaw, 8 1/2" (21.5cm)



PG357622 : YASARGIL VASCULAR SCISSORS

Straight, 7" (18.0cm), 12.0mm blades





PG357623 : JACOBSON VASCULAR SCISSORS

Curved, 7" (18.0cm), 12.0mm blades



PG357642 : YOU-POTTS VASCULAR SCISSORS

Narrow, thin blade and round spring handles designed for a more precise and controlled arteriotomy. Angled 90, 7" (18.0 cm), 14.0 mm blades



PG357643 : YOU-POTTS VASCULAR SCISSORS

Narrow, thin blade and round spring handles designed for a more precise and controlled arteriotomy. Angled 120, 7" (18.0 cm), 12.0 mm blades



PG357645 : MICRO RING TIP FORCEPS

Smooth, 1.5mm tip, 7" (18 cm)



PG357645DG : MICRO DIAMA-GRIP RING FORCEPS

1.5 mm tip, 7" (18 cm), Diama-Grip jaws

