




PG359307 : CRILE TITANIUM HEMOSTATIC FORCEPS

Curved, Serrated, 5 5/8" (15 cm)



PG359313 : O'SHAUGHNESSY TITANIUM FORCEPS

Right angle, Serrated, 8 1/4" (21 cm)




PG359314 : HEISS TITANIUM FORCEPS

Angled, Serrated, 8 1/4" (21 cm)




PG359321 : O'SHAUGHNESSY TITANIUM FORCEPS

Right angle, Serrated, 2.5mm tip, 9 1/2" (24 cm)




PG359325 : CRAFOORD TITANIUM CLAMP FORCEPS

Curved, Serrated, 9 1/2" (24 cm)



PG359328 : CASTROVIEJO MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, Without lock, 1.1mm Diama-Grip jaw, 9" (23cm)



PG359329 : CASTROVIEJO MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 1.1mm Diamo-Grip jaw, 9" (23cm)




PG359335 : MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 9 x 1.2mm jaw, 7 1/8" (18cm), TC



PG359339 : MICRO TITANIUM NEEDLE HOLDER

Straight, Flat handle, With lock, 1.4mm jaw, 9" (23cm), TC



PG359340 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, Without lock, 1.1mm Diamo-Grip jaw, 9" (23cm)



PG359342 : MICRO TITANIUM NEEDLE HOLDER

Curved, Round handle, Without lock, 1.1mm Diamo-Grip jaw, 9" (23cm)



PG359343 : MICRO TITANIUM NEEDLE HOLDER

Curved, Round handle, With lock, 1.1mm Diamo-Grip jaw, 9" (23cm)



PG359344 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, Without lock, 1.1mm jaw, 6 1/4" (16cm), TC

PG359345 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, With lock, 1.1mm jaw, 6" (16cm), TC

NO IMAGE
AVAILABLE


PG359346 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, Without lock, 7 1/8" (18cm), TC




PG359347 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, With lock, 8.8 x 1.3mm jaw, 7 1/8" (18cm), TC




PG359349 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, Without lock, 1.5mm jaw, 8 1/4" (21cm), TC



PG359350 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, With lock, 1.4mm jaw, 8 1/4" (21cm), TC



PG359352 : MICRO TITANIUM NEEDLE HOLDER

Straight, Round handle, With lock, 1.5mm jaw, 9" (23cm), TC



PG359353 : PEDIATRIC TITANIUM STERNAL SPREADER


Spread 2" (5.0cm), Arm length 52mm, 3" (7.5cm)






PG359354 : DEBAKEY TITNIUM PEDIATRIC CLAMP

Straight, 1.6mm tip width, 4 3/8" (11cm)




PG359355 : DEBAKEY TITNIUM PEDIATRIC CLAMP

Angled 60, 7/8" (2.3cm) tip length, 1.6mm tip width, Curved shaft, 4" (10.5cm)



PG359356 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Angled, 3/8" (0.95cm) tip length, 1.6mm tip width, Angled shaft, 3 7/8" (10cm)



PG359357 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Spoon, 1/4" (0.64cm) tip length, 1.2mm tip width, 3 7/8" (9.5cm)



PG359358 : DEBAKEY TITNIUM ANASTOMOSIS CLAMP

Angled, 5/8" (1.59cm) tip length, 1.6mm tip width, 4 1/8" (10.5cm)

