

GENERAL INSTRUMENTS



PO271020 : HEPARIN CANNULA FINE MODEL



PO271021 : SUCTION TUBE NOWACK 23CM DIAM 4MM



PO271052 : SUCTION TUBES CONE-LUER ANGLED D=0,9MM



PO271053 : SUCTION TUBE CONE LUER 10CM ANGLED DIAM 1.2MM



PO271055 : SUCTION TUBES CONE-LUER ANGLED L=11,5CM D=1,8MM



PO271056 : SUCTION TUBES CONE-LUER ANGLED L=10CM D=1,9MM



PO271057 : SUCTION TUBE CONE LUER 10CM ANGLED DIAM 2.4MM



PO271067 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM D=0,9MM



PO271068 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM D=1,5MM



PO271069 : SUCTION TUBES CURVED OLIVE SHAPED TIP L=12CM D=1,9MM



PO271070 : SUCTION TUBE 12CM CURVED OLIVE SHAPED TIP 2.5MM



PO271071 : SUCTION TUBE 12CM CURVED OLIVE SHAPED TIP 3MM



PO271072 : SUCTION TUBE 12CM CURVED OLIVE SHAPED TIP 4MM

PO271089 : SUCTION TUBES VEIN CANNULA L=8CM D=1,2MM



PO271090 : SUCTION TUBES VEIN CANNULA L=9,5CM D=1,5MM



PO273001 : TROCAR PLEUROTOME BLUNT TIP 18



PO273003 : TROCAR PLEUROTOME BLUNT TIP 22



PO273004 : TROCAR PLEUROTOME BLUNT TIP 24




PO273005 : TROCAR PLEUROTOME POINTED TIP 16



PO273006 : TROCAR PLEUROTOME POINTED TIP 18





PO273007 : TROCAR PLEUROTOME POINTED TIP 20



PO273008 : TROCAR PLEUROTOME POINTED TIP 22



PO273009 : TROCAR PLEUROTOME POINTED TIP 24



PO281000 : SUTURE INST NEEDLES REVERDIN INTESTINAL STRAIGHT 15CM



PO281001 : SUTURE INST NEEDLES REVERDIN INTESTINAL HALF CURVED 15CM

