

GENERAL**PG480105 : INFANT RETRACTOR**

1/2" (13.5mm) width x 1/2" (1.0cm) depth blade, 5 5/8" (14.0cm)

**PG480140 : ROUX RETRACTOR, SMALL**


Double ended, 1 1/2" (3.5cm) / 1 1/8" (2.5cm) width blade, 7" (18.0cm)

**PG480145 : ROUX RETRACTOR, LARGE**

Double ended, 1 3/4" (4.5cm) / 1 3/8" (3.5cm) width blade, 7 1/8" (18.0cm)

**PG480170 : PARKER RETRACTOR, SMALL**

Double ended, 5/8" (1.5cm) / 5/8" (1.5cm) width blade, 5 3/8" (14.0cm), 2/set

**PG480172 : PARKER RETRACTOR, LARGE**

Double ended, 1" (2.5cm) / 1" (2.5cm) width blade, 7 1/8" (18.0cm), 2/set

**PG480195 : RICHARDSON EASTMAN RETRACTOR**

Set of 2. Large 10 5/8" (27cm), 63mm x 42mm and 39mm x 37mm blades, Small 9 7/8" (20cm), 30mm x 29mm and 24mm x 22mm blades.




PG480196 : RICHARDSON EASTMAN SMALL RETRACTOR

7/8" (2.0cm) / 1 1/8" (3.0cm) width blade, 9 7/8" (25.0cm)



PG481220 : DECOMPRESSION CUSHING RETRACTOR

10mm width blade, 8 1/8" (20.5cm)



PG481231 : CUSHING VAIN / NERVE RETRACTOR

5/8" (14.5mm) width x 1/2" (1.5cm) depth blade, 12 3/4" (32.5cm)



PG481500 : LITTLE RETRACTOR

1/2" (12mm) width x 3/4" (2.0cm) depth blade, 7 1/2" (19.0cm)



PG481520 : GREEN RETRACTOR

21mm width x 1" (2.5cm) depth blade, 8 7/8" (22.5cm)



PG481530 : MCBURNEY RETRACTOR

5/8" (16mm) width x 3/4" (2.0cm) depth blade, 7 5/8" (19.5cm)



PG481610 : LAHEY RETRACTOR

5mm width x 1 1/8" (3.0cm) depth blade, 7 3/4" (19.5cm)

PG481617 : RICHARDSON RETRACTOR

Hollow handle, 1.25" (27mm) width x 1 1/4" depth blade, 9 1/2" (24.0cm)



PG481618 : RICHARDSON RETRACTOR

Hollow handle, 1.5" (37mm) width x 1 1/2" depth blade, 9 1/2" (24.0cm)



PG481629 : RICHARDSON RETRACTOR

Hollow handle, 0.75" (21mm) width x 2" depth blade, 9 1/2" (24.0cm)



PG482371 : HARRINGTON RETRACTOR

Hollow handle, 65mm width blade, 13" (33.0cm)




PG488306 : WEITLANER RETRACTOR

2 x 3 Blunt prongs, 4 1/2" (11.5cm)



PG488361 : ADSON RETRACTOR

Spread 90mm, 4 x 4 Blunt prongs, 7 1/2" (19cm)



PG496105 : LAHEY DISSECTING FORCEPS

Strongly Curved, 22 x 2.4mm tip, 7 7/8" (20 cm)






PG496111 : LAHEY DISSECTING FORCEPS

Strongly Curved, 2.5mm tip, 7 3/4" (20 cm)




PG496216 : STOREY HEMOSTAT FORCEPS

Strongly Curved, Serrated, 8 5/8" (22 cm)



PG496217 : STOREY HEMOSTATIC FORCEPS

Right angle, Serrated, 8 3/8" (22 cm)



PG536220 : ANDREWS PYNCHON SUCTION TUBE

Curved, 6mm tip, 3mm dia. removable shaft, 9 5/8" (24.5cm)



PG640280 : POTTS TENOTOMY SCISSORS

Curved, 6" (15.0cm)

