


GENERAL




PG166495 : WILLIAMS LACRIMAL PROBE

Size 7-8, 5 1/8" (13cm)




PG166500 : PROBES WITH PROXIMAL EYE

Stainless Steel, 4 1/2" (11.5cm)




PG166501 : PROBES WITH PROXIMAL EYE

Stainless Steel, 5 1/8" (13cm)



PG166502 : PROBES WITH PROXIMAL EYE

Stainless Steel, 5 3/4" (14.5cm)




PG166503 : PROBES WITH PROXIMAL EYE

Stainless Steel, 5 7/8" (15cm)




PG166504 : PROBES WITH PROXIMAL EYE

Stainless Steel, 7 1/8" (18cm)



PG166505 : PROBES WITH PROXIMAL EYE

Stainless Steel, 7 7/8" (20cm)



PG166507 : PROBES WITH PROXIMAL EYE

Stainless Steel, 9 7/8" (25cm)



PG166633 : GROOVED DIRECTOR, PROBE POINT

Stainless Steel, 6 1/4" (16cm)



PG181220 : ADSON DRESSING FORCEPS

Serrated, 4 3/4" (12 cm)



PG181223 : ADSON TISSUE FORCEPS

1 x 2 Teeth, 4 3/4" (12 cm)



PG181225 : ADSON TISSUE FORCEPS

2 x 3 Teeth, 5" (13 cm)



PG181228 : ADSON DRESSING FORCEPS

Serrated, 1 x 2 Teeth, 4 3/4" (12 cm)

PG181235 : BROWN-ADSON TISSUE FORCEPS

9 x 9 Teeth, 4 3/4" (12 cm)



PG181237 : EWALD DRESSING FORCEPS

Serrated, 4 3/4" (12 cm)



PG181240 : EWALD TISSUE FORCEPS

1 x 2 Teeth, 4 3/4" (12 cm)



PG181255 : BROWN TISSUE FORCEPS

Straight, 9 x 9 Teeth, 6" (16 cm)



PG181360 : RUSSIAN TISSUE FORCEPS

6" (15 cm)



PG181365 : RUSSIAN TISSUE FORCEPS

8" (21 cm)



PG181370 : RUSSIAN TISSUE FORCEPS

10" (26 cm)



PG181410 : CUSHING DRESSING FORCEPS

Serrated, 7" (18 cm)



PG183010 : BABCOCK TISSUE FORCEPS

6 1/4" (16 cm)



PG183020 : BABCOCK TISSUE FORCEPS

7 7/8" (20 cm)



PG342965 : BROWN TISSUE FORCEPS

9 x 9 Teeth, 8" (21 cm)



PG343395 : THOMPSON THORACIC INSTRUMENT HOLDER

