


# TAUT® SYSTEM ONE™

## Comprehensive laparoscopic common bile duct exploration kit

The System ONE is an integrated combination of products designed specifically for t[illegible]ystic laparoscopic com[illegible] bile duct exploration. The tools needed for several methods to clear CBD stones are provided in one convenient kit.

### MULTIPURPOSE BALLOON CATHETER

- **FLUSH** Flush small stones and debris through central lumen
- **EXPLORE** Insert Titanium Retrieval Basket through central lumen for fluoroscopic exploration
- **DILATE** Control dilation and ease guidewire placement for choledochoscopy

### TITANIUM4™ RETRIEVAL BASKET

- Fully articulated four-wire titanium alloy basket can be rotated 360°
- Separate controls for deployment and rotation on ergonomic handle
- Resists breaking and deforming

Teleflex MEDICAL

## SYSTEM ONE™ COMPREHENSIVE LAP CBDE KIT

TAUT®

| REF NUMBER | DESCRIPTION |
|---|---|
| 50000 | System ONE Reusable Organizer |
| | Titanium⁴ Retrieval Basket |
| | Multipurpose Balloon Catheter |
| | Percutaneous Intraducer |
| | Softwire Guidewire |
| | High Pressure Syringe |




SYSTEM ONE KIT

## DISPOSABLE STERILE COMPONENTS

**All items can be ordered to restock kit**

### MULTIPURPOSE BALLOON CATHETER (1 UNIT PER BOX)

TAUT®

| REF NUMBER | DESCRIPTION | CATHETER DIAMETER | INFLATED DIAMETER | LENGTH |
|---|---|---|---|---|
| 50640 | Balloon catheter with radiopaque markers for clear visualization of distal and proximal balloon ends | 5 French/1.67mm | 6 mm | Balloon 4 cm<br>Total length 75cm |



MULTIPURPOSE BALLOON CATHETER

### TITANIUM4™ RETRIEVAL BASKET (1 UNIT PER BOX)

TAUT®

| REF NUMBER | DESCRIPTION | WIRE DIAMETER | OPEN BASKET DIAMETER | LENGTH |
|---|---|---|---|---|
| 50024 | 4 Flat wire titanium alloy basket with 360 degree articulation | 2.4 French/0.8mm | 15mm | 100cm |



TITANIUM⁴ RETRIEVAL BASKET

### PERCUTANEOUS INTRADUCER® (10 UNITS PER BOX)

TAUT®

| REF NUMBER | DESCRIPTION | DIAMETER | LENGTH |
|---|---|---|---|
| PI-128 | Accomodates wire basket, balloon catheter, guidewire, and many choledochoscopes | 10 French/3.33mm | 20.3cm/8in |

### SOFTWIRE™ GUIDEWIRE (1 UNIT PER BOX)

TAUT®

| REF NUMBER | DESCRIPTION | DIAMETER | LENGTH |
|---|---|---|---|
| 50035 | Dual flexible ended PTFE Teflon® coated guidewire | 0.035in | 125cm/49.2in |

### HIGH PRESSURE SYRINGE (1 UNIT PER BOX)

TAUT®

| REF NUMBER | DESCRIPTION |
|---|---|
| K05-01251 | High pressure syringe for inflation of balloon catheter |



HIGH PRESSURE SYRINGE

**TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990**

**TELEFLEX MEDICAL** Research Triangle Park, NC 27709
Phone: 866.246.6990 Fax: 800.399.1028 Intl: 919.361.3961
**WWW.TELEFLEXMEDICAL.COM**
Taut, Softwire, Intraducer, Titanium⁴ and System ONE are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. Teflon is a registered trademark of E.I. du Pont de Nemours and Company.
©2010 Teleflex Incorporated. All rights reserved. 2009-0214 Rev 00

Teleflex MEDICAL