WECK®



Confidence
Clarity
Control

Teleflex®

FROM WECK®. THE MARKET LEADER IN LIGATION, COMES A NEW LINE FOR ACCESS.



WECK Vista™
BLADELESS LAPAROSCOPIC ACCESS PORTS

OPTICAL
- Clear visualization of tissue layers
- Securely grasps laparoscopes

BALLOON/CONE
- Designed to maintain pneumoperitoneum
- Easy one-hand stabilizer mechanism

BLADELESS PORTS
- Asymmetrical dilating tip
- Ridged or smooth cannulas

BALLOON/CONE **BLADELESS PORTS**





REUSABLE OBTURATORS
- Bladeless tip compatible with cannula only offerings
- Easy to reprocess
- Cost-effective solution

ORDERING INFORMATION

WECK Vista PORTS, 5 MM

RIDGED CANNULAS REF. NO.	SMOOTH CANNULAS REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
43507R	43507	WECK Vista Port	3 – 5mm X 70mm	6
40507R	40507	WECK Vista Port	5mm X 70mm	6
40511R	40511 (Not Optical)	WECK Vista Optical Port	5mm X 100mm	6
40511RC	40511C	WECK Vista Cannula	5mm X 100mm	6
40513R	40513	WECK Vista Port	5mm X 125mm	6

WECK Vista PORTS, 8 MM

RIDGED CANNULAS REF. NO.	SMOOTH CANNULAS REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
40811R	N/A	WECK Vista Port	5 – 8mm X 100mm	6

WECK Vista PORTS,10 MM

RIDGED CANNULAS REF. NO.	SMOOTH CANNULAS REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
405910R	405910 (Not Optical)	WECK Vista Optical Port	5 – 10mm X 100mm	6
405910RC	405910C	WECK Vista Cannula	5 – 10mm X 100mm	6

WECK Vista PORTS, 12 MM

RIDGED CANNULAS REF. NO.	SMOOTH CANNULAS REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
405912R	405912 (Not Optical)	WECK Vista Optical Port	5 – 10 – 12mm X 100mm	6
405912RC	405912C	WECK Vista Cannula	5 – 10 – 12mm X 100mm	6
40591213R	40591213	WECK Vista Port	5 – 10 – 12mm X 125mm	6

WECK Vista PORTS, 15 MM

RIDGED CANNULAS REF. NO.	SMOOTH CANNULAS REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
41510R	41510	WECK Vista Port	15mm X 110mm	6

WECK Vista BALLOON ACCESS PORTS

REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
410944S	WECK Vista Balloon Access Port	5 – 10mm X 53mm	6
410944	WECK Vista Balloon Access Port	5 – 10mm X 70mm	6
410944L	WECK Vista Balloon Access Port	5 – 10mm X 100mm	6
412944	WECK Vista Balloon Access Port	5 – 10 – 12mm X 70mm	6
412944L	WECK Vista Balloon Access Port	5 – 10 – 12mm X 100mm	6

WECK Vista CONE OPEN ACCESS PORTS

REF. NO.	DESCRIPTION	SIZE	QTY PER BOX
405933	WECK Vista Cone Open Access Port	5 – 10 – 12mm X 100mm	6
405933L	WECK Vista Cone Open Access Port	5 – 10 – 12mm X 125mm	6

INSUFFLATION NEEDLES

REF. NO.	DESCRIPTION	LENGTH	QTY PER BOX
35120	Insufflation Needle 14 Gauge	120mm	10
35150	Insufflation Needle 14 Gauge	150mm	10

WECK REUSABLE OBTURATORS

REF. NO.	DESCRIPTION	SIZE	COMPATIBLE SINGLE USE CANNULAR	QTY PER BOX
4050788S	WECK Reusable Obturator	5mm X 70mm	4050722	1
405188	WECK Reusable Obturator	5mm X 100mm	40511C, 40511RC	1
45188S	WECK Reusable Obturator	5mm X 100mm	451022U	1
410988	WECK Reusable Obturator	10mm X 100mm	405910C, 405910RC	1
412988	WECK Reusable Obturator	12mm X 100mm	405912C, 405912RC	1

For a complete description of all indications and contraindications, please refer to the associated Instructions for Use.

http://teleflex.orlive.com
www.weckvista.com

Teleflex, Weck and WECK Vista are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

© 2012 Teleflex Incorporated. All rights reserved. 2012-1206

Teleflex®

TELEFLEX MEDICAL PO Box 12600 Research Triangle Park, NC 27709
Toll Free: 866.246.6990 Phone: 919.544.8000 Intl: 919.433.8088
TELEFLEX.COM