WECK®

WECK® REUSABLE OBTURATORS

Confidence. Clarity. Control.



- Controlled entry in a durable stainless steel design
- Bladeless asymmetrical dilating tip
- Easy to reprocess one piece design

Teleflex®

WECK® REUSABLE OBTURATORS



The Weck Reusable Obturators offer convenience and are eco-friendly which reduce waste in the OR. The easy to reprocess one piece design is sturdy to offer surgeons a controlled entry while the bladeless asymmetrical tip spreads tissue in a natural fashion.

WECK Vista™ single-use cannulas are compatible with the reusable obturators. Conical entry for easy instrument exchange is offered on some models as well as ridged or smooth cannulas.



BENEFITS OF REUSABLE OBTURATORS:

1. Asymmetrical tip designed to minimize tissue trauma
2. Eco-friendly which reduces waste in the OR
3. Low maintenance one piece design
4. Sturdy stainless steel for controlled entry

ORDERING INFORMATION

WECK REUSABLE OBTURATORS

REF. NO.	DESCRIPTION	SIZE	COMPATIBLE SINGLE-USE CANNULA	QTY PER BOX
4050788S	WECK Reusable Obturator	5 mm X 70 mm	4050722	1
405188	WECK Reusable Obturator	5 mm X 100 mm	40511C, 40511RC	1
45188S	WECK Reusable Obturator	5 mm X 100 mm	451022U	1
410988	WECK Reusable Obturator	10 mm X 100 mm	405910C, 405910RC	1
412988	WECK Reusable Obturator	12 mm X 100 mm	405912C, 405912RC	1

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Teleflex, Weck and WECK Vista are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
©2012 Teleflex Incorporated. All rights reserved. 2012-1372

Teleflex®

TELEFLEX MEDICAL PO Box 12600 Research Triangle Park, NC 27709
Toll Free: 866.246.6990 Phone: 919.544.8000 Intl: 919.433.8088
TELEFLEX.COM