

Pilling®



Improved Tucker



Carlens



Pilling®
Mediastinoscopy
Basic Instrument Configuration

Designed for Comprehensive Mediastinal Visualization


- Pilling® Carlens and Tucker Mediastinoscopes available in five sizes
- Select the mediastinoscopes, instrumentation, and tray configuration to best suit your needs for procedural optimization

For more information, please contact your Pilling® Sales Representative or visit our product catalog at
www.teleflexsurgicalcatalog.com


Teleflex®

Drawings illustrate Nominal Proximal Internal Dimensions and Minimum Internal Diameter (I.D.), **NOT TO SCALE**

Carlens



Improved Tucker



Mediastinoscopes

All scopes supplied with light carrier.

Carlens

☐	522300P	Short, 16 mm x 116 mm
☐	522301	Long, 18 mm x 162 mm
☐	522302	Extra-Long, 19 mm x 191 mm

Improved Tucker

☐	522306	Short, 24 mm x 147 mm
☐	522311	Long, 27 mm x 180 mm

Instruments and Accessories

Forceps

☐	505100	Jackson 4 mm Straight Cup Forceps, 28 cm
☐	505104	Jackson 4 mm Angled Up Cup Forceps, 28 cm
☐	505110	Jackson 2 mm Straight Cup Forceps, 20 cm
☐	505112	Jackson 2 mm Straight Cup Forceps, 28 cm
☐	505114	Jackson 2 mm Angled Cup Forceps, 20 cm
☐	505116	Jackson 2 mm Angled Cup Forceps, 28 cm
☐	505230	Patterson 4 mm Elongated Cup Forceps, 28 cm
☐	505400	Jackson Grasping Forceps, 28 cm
☐	505420	Jackson Delicate Grasping Forceps, 30 cm
☐	506747	Norris Sponge Holding Forceps
☐	506748	Mediastinal Dissecting Forceps

Aspiration

☐	507453	Aspirating Needle, 23 Gauge, 14.7 mm, 21 cm Working Length
☐	507454	Spare Distal Needle Tip

Suction

☐	507045	Insulated Suction Tube, 12 fr, 22 cm, with 10 ft Bovie Cord
☐	351970	Bovie Electrosurgical Handle, 13 cm

Instruments and Accessories

Ligation

☐	506743	McKenzie Sterling Silver "V" Ligating Clips, Smooth, 100/Vial
☐	506744	McKenzie Sterling Silver "V" Ligating Clips, Serrated, 100/Vial
☐	506746	Duane Sterling Silver "U" Ligating Clips, Smooth, 100/Vial
☐	506737	Drew Alligator Clip Applying Forceps, Serrated
☐	506742	Scoville Alligator Clip Applying Forceps, Slotted
☐	506749	Clip Holding Magazine

Illumination

A Blue Fiber Optic Cable is required for connecting Carlens Mediastinoscopes to lightsource.

☐	528305XLTP	Xenalight® XLT II Lightsource, 5-Port Turret
☐	P358PRA	Pilling® to Pilling® Blue Fiber Optic Cable, 8 ft
☐	A358PRA	ACMI to Pilling® Blue FIber Optic Cable, 8 ft
☐	50616AC	Pilling® to ACMI Illuminator-End Adapter for Tucker Light Cable

Holders

☐	502248	Tobey Laryngoscope Holder (Compatible with Improved Tucker Mediastinoscope)

Tray Setup

☐	P37331	Single-Tier ENT Tray w/o Brackets, 11" x 23" x 3.5"
☐	P37338	Tray Shelf, 10.5" x 21.5" x 2.25"
☐	P37545 x4	"L" Brackets to Support Shelf
☐	P37306 x3	Universal Keyhole Bracket 15 Position, 10"
☐	P37379 x5	Keyhole Bracket, Single, 1.5"
☐	087396 x2	Ultrasoft Silicone Spiked Mat, 11" W x 22" L

Additional Needs

Recommended Configuration(s) Selected and Approved by:

Print Name **Date**

Signature

Sales Rep Name **Sales Rep Phone Number**

Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are property of their respective owners. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 02/2020. © 2020 Teleflex Incorporated. All rights reserved. MC-006131

Teleflex®

Teleflex · 3015 Carrington Mill Boulevard, Morrisville, NC 27560
Toll Free: 866 246 6990 | Phone: +1 919 544 8000
teleflex.com