
Pilling

Extended reach
for a range of patient sizes and anatomies



Light carrier
is flush and parallel to the handle to help keep light cables out of the way

Scope channel shape
helps visualize surgical space and allows instrument pass-through

New Pilling Laryngoscope!

Pilling®

Garrett-Ossoff-Pilling
Laryngoscope

Garrett-Ossoff-Pilling Laryngoscope

The new Garrett-Ossoff-Pilling Laryngoscope helps surgeons focus on critical anatomy during procedures in a range of patient sizes, and is designed for patients with a greater distance between their teeth and vocal cords.

Teleflex®

The Ossoff-Pilling Laryngoscope Family:




Laser-compatible scope allows you to have binocular vision.

522196
Garrett-Ossoff-Pilling
Adult - Extra-long

522190
Ossoff-Pilling
Adult Female

522191
Ossoff-Pilling
Adult Male

Additional Pilling Laryngoscope Designs:



Ossoff-Pilling Laryngoscopes

CATALOG #	SIZE	OVERALL LENGTH	MINIMAL DISTAL OPENING	MOUTH PROXIMAL DIMENSIONS	REPLACE CARRIER
522196	Extra Long	188 mm	8.0 mm	25.0 mm x 15.0 mm	527868
522191	Male	168 mm	8.0 mm	25.0 mm x 15.0 mm	527865
522190	Female	168 mm	7.0 mm	19.0 mm x 15.0 mm	527864

teleflexsurgicalcatalog.com

Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (12/2016). Subject to technical changes without further notice. © 2016 Teleflex Incorporated. All rights reserved. MC-002813

Teleflex®

Teleflex · 3015 Carrington Mill Boulevard, Morrisville, NC 27560
Toll Free: 866 246 6990|Phone: +1 919 544 8000
teleflex.com