

TAUT® CHOLANGIOGRAM CATHETERS & INTRADUCERS®

Teleflex MEDICAL


THE TAUT ADVANTAGE

METAL SUPPORT TUBE
Maximize clip pressure without occluding the catheter to create a tight seal and prevent extravasation

MEMORY
Customizable tip angle for efficiency and ease of insertion into the cystic duct

MULTIPLE LENGTHS
Available in 18 in. (45.7 cm) and 30 in. (76 cm)

CHOLANGIOGRAM CATHETER INTRADUCERS
Provide access to the cystic duct for cholangiography and laparoscopic common bile duct exploration



Catheter memory eases insertion

CHOLANGIOGRAM CATHETER TIP DESIGNS



20018-M55



20018-M56



20018-M57



20018-010



18305



18400

CHOLANGIOGRAM CATHETERS

All products packed sterile, 10 units/box.

Product Number	Catheter Diameter	Tip Diameter	Length	Metal Support Tube
20018-M55	4.5 French/1.5 mm	4.5 French/1.5 mm	18in./45.7 cm	Yes
20018-M56	4.5 French/1.5 mm	5.5 French/1.8 mm	18in./45.7 cm	Yes
20018-M57	4.5 French/1.5 mm	6.5 French/2.2 mm	18in./45.7 cm	Yes
20018-010	4.5 French/1.5 mm	7.5 French/2.5 mm	18in./45.7 cm	No
18305	4.5 French/1.5 mm	4.5 French/1.5 mm	30in./76.0 cm	No
18400	4.5 French/1.5 mm	4.5 French/1.5 mm	30in./76.0 cm	Yes

INTRADUCERS

All products packed sterile, 10 units/box.

Product Number	Internal Diameter	Length	Compatible Catheters
PI-63	4.5 French/1.6 mm	3.0in./7.6 cm	20018-M55, 18305, 18400
PI-93*	7.5 French/2.4 mm	3.5in./8.9 cm	20018-M55, 20018-M56, 20018-M57, 18305, 18400
PI-104	8.0 French/2.7 mm	3.5in./8.9 cm	20018-M55, 20018-M56, 20018-M57, 20018-010, 18305, 18400
PI-128	10 French/3.3 mm	8.0in./20.3 cm	20018-M55, 20018-M56, 20018-M57, 20018-010, 18305, 18400

* Contains natural latex rubber. All other items are latex-free.

For more information on Taut products, please call contact your local Teleflex Medical surgical specialist, or call customer service at 866.246.6990.

Taut and Intraducer are registered trademarks of Teleflex Incorporated or its affiliates.
©2009 Teleflex Incorporated Lit. No. 2009-0075 Rev. 01 Research Triangle Park, NC
Tel: 866.246.6990 Intl: 919.544.8000 www.teleflexmedical.com

Teleflex MEDICAL